

```

1  REM echo off
  @REM /* -----
  @REM Covidien Copyright 2008
  @REM
5  @REM File Name: 10039032 PB560 Pediatric Volume and Pressure Test Script
  @REM
  @REM Function: Script file to perform PB560 Pediatric Patient Circuit Tes
  @REM This file is executed by the OpenVent s/w tool p/n 10026042
  @REM
10 @REM
  @REM History:
  @REM
  @REM Revision Date Author Comment
  @REM X1 12-June-09 B.Hurd Initial Version from PB540 10025468 Re
15 @REM A 23-Feb-10 A.Donnelly Update to rev A for P4
  @REM No functional change to tests
  @REM
  @REM B 15-March-12 B.Hurd Add Targeted Vent for OpenVent 1.6.0.7
  @REM universal test PC and Map T:\ to
20 @REM \\galwaylapo3\TESTRES$\PB560\%COMPUTERNAME%
  @REM
  @REM C 30-Nov-18 B.Hurd/ C. Chiranjil Changed I/E to Insp Time for all the c
  @REM with CPU S/W. Correct test titles for Insp Ti
  @REM to be in seconds.
25 @REM Update to new Insp Sens settings: OFF1=0p,2=1p,3=2
  @REM 4=3,5=4,6=5.
  @REM Changed %TI_SETTING% to I_ADJUST_ADJUST_TI_CONTROL
  @REM Update alarm compare commands to use parms
  @REM HIGH_PRESSURE_ALARM and DISCONNECTION_ALARM
30 @REM for test report readability. These parms default
  @REM to zero.
  @REM Update test limits for monitored PEEP and monitored Pi
  @REM from +/- (2+8%) to +/- (2+4%) to align with new Homecare
  @REM According to EN 80601-2-72:2015 Clause 201.12.4.102,
35 @REM Indicated Airway pressure accuracy shall be within ± (
  @REM Update SafetyPressure from 60 to 90 to match change pe
  @REM SCR195 CDP-00014971 which has been made to align with
  @REM Homecare Standard (EN 80601-2-72:2015) Clause 201.12.4
  @REM -----
40 @REM
  @REM ----- References -----
  @REM [1] PB540/PB560 Volume and Pressure Tests, Manuf Test Spec Doc, P/N 10025
  @REM [2] OpenVent Software Specification, P/N 10026043
45 @REM ----- Comparison Operators -----
  @REM
  @REM EQU - Equal To
  @REM NEQ - Not Equal To
50 @REM LSS - Less Than
  @REM LEQ - Less Than or Equal To
  @REM GTR - Greater Than
  @REM GEQ - Greater Than or Equal To
  @REM
55 @REM -----
  @REM -----
  @REM :CHK_DRV_MAP560
60
  ECHO CHECK FOR MAPPED NETWORK DRIVE T:\PB560
  IF EXIST T:\PB560 GOTO CHK_DRV_MAP560_END
65 @REM delete previous T:\ mapping to other PB model
  net use T: /DELETE
  ECHO T:\DRIVE IS NOT MAPED
  ECHO MAPPING DRIVE T:\ to \\galwaylapo3\TESTRES$\PB560\%COMPUTERNAME%
70
  net use T: \\galwaylapo3\TESTRES$\PB560\%COMPUTERNAME%

```

```

      :CHK_DRV_MAP560_END
@REM -----
75  @REM ----- Constants -----
    set output_dir="T:\PB560\Final Test\Ped"
    set report_name="PB560 Pediatric Volume and Pressure Tests"
    set compliance=C3
80  set resistance=R50

    @REM ----- Test Program Link -----
    @REM This line inserts the link to program testventil.exe
    @REM %TestProgram% represents the path for testventil.exe
85  set TestProgram="C:\Program Files\OpenVent\testventil\testventil.exe"
    set GenerateProgram="C:\Program Files\Openvent\GenerateReport.exe"

    REM ----- Targeted vent -----
    SET pTargetedVent=40966
90

    @REM ----- Creation of the XML output file -----
    REM add 'draw' at end of command to generate curves for report
    call %TestProgram% -F %output_dir% %report_name% %compliance% %resistance%
    IF %ERRORLEVEL% NEQ 0 GOTO end
95

    REM ----- Zero PTS2000 Low Pres Input -----
    call %TestProgram% -Z
    IF %ERRORLEVEL% NEQ 0 GOTO end
100

    REM ----- Ped Circuit 'Yes' -----
    call %TestProgram% -R CONFIG_ADJUST_PEDIATRIC_CIRCUIT 1
    IF %ERRORLEVEL% NEQ 0 GOTO end
105

    REM ----- Relative pressure 'Yes' -----
    call %TestProgram% -R CONFIG_PRESSURE_SUPPORT_RELATIVE 1
    IF %ERRORLEVEL% NEQ 0 GOTO end

110  REM ----- VOL AC 50mL Test -----

    SET title="Mode= VOL, Vol=50mL, Rate=30, Insp time= 0.7s, Peep=0, Insp Sen

115  REM modes: VOL=0 PRES=1 PSVT=2 VSIMV=3 CPAP=4
    call %TestProgram% -R CONFIG_ADJUST_MODE 0
    IF %ERRORLEVEL% NEQ 0 GOTO end

120  REM Set PIP widest alarm limits 20-600 (note: 20 equals 2)
    call %TestProgram% -R ADJUST_ADJUST_LOW_PIP 20
    IF %ERRORLEVEL% NEQ 0 GOTO end
    call %TestProgram% -R ADJUST_ADJUST_HIGH_PIP 600
    IF %ERRORLEVEL% NEQ 0 GOTO end
125

    REM ----- Set Alarm Limits to OFF -----
    call %TestProgram% -R ADJUST_LOW_VTE_NO_SELECT 1
    IF %ERRORLEVEL% NEQ 0 GOTO end
    call %TestProgram% -R ADJUST_HIGH_VTE_NO_SELECT 1
130  IF %ERRORLEVEL% NEQ 0 GOTO end
    call %TestProgram% -R CONFIG_FIO2_HIGH_SELECT 0
    IF %ERRORLEVEL% NEQ 0 GOTO end
    call %TestProgram% -R CONFIG_FIO2_LOW_SELECT 0
    IF %ERRORLEVEL% NEQ 0 GOTO end
135

    REM Turn OFF Barometric Compensation (0 = OFF)
    call %TestProgram% -R CONFIG_ADJUST_BAROMETRIC_CORRECTOR 0
    IF %ERRORLEVEL% NEQ 0 GOTO end

140  REM set PEEP to 0
    call %TestProgram% -R ADJUST_PEEP_NO_SELECT 0
    IF %ERRORLEVEL% NEQ 0 GOTO end
    call %TestProgram% -R ADJUST_ADJUST_PEEP 0
    IF %ERRORLEVEL% NEQ 0 GOTO end

```

```

145      REM RAMP 0 = Slope, RAMP 1 = Rectangle
      call %TestProgram% -R ADJUST_ADJUST_RAMP 0
      IF %ERRORLEVEL% NEQ 0 GOTO end

150      REM Insp Time
      call %TestProgram% -R ADJUST_ADJUST_TI_CONTROL 700

      REM allow insp sens trigger setting to 5
155      call %TestProgram% -R ADJUST_ADJUST_INSP_SENS 6
      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM Set R-Rate
      call %TestProgram% -R ADJUST_ADJUST_CONTROL_R 30
160      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM Set volume (mL)
      call %TestProgram% -R ADJUST_ADJUST_VOL_CONTROL 50
      IF %ERRORLEVEL% NEQ 0 GOTO end

165

      REM ----- REFRESH VENTILATOR SCREEN -----
      REM (cursor needs to be on the first setting at the top of the main menu i

170      call %TestProgram% -R KEYBOARD_KEYBOARD_UP_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
175      REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_DOWN_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
180      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
185      call %TestProgram% -R KEYBOARD_KEYBOARD_UP_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
190      REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_DOWN_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM -----

195      REM --- SETTINGS COMPARISONS--- compares Vent settings to PTS2000 measurem

      SET peep_tolerance=1.5
      SET VOL_SETTING=50
      REM Vti/Vte Tolerance is 10% + 10ml

200      call %TestProgram% -C %VOL_SETTING% VTI 10p+10
      IF %ERRORLEVEL% NEQ 0 GOTO end
      call %TestProgram% -C I_ADJUST_ADJUST_CONTROL_R RATE 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
205      call %TestProgram% -C I_ADJUST_ADJUST_PEEP/10 PEEP %peep_tolerance%
      IF %ERRORLEVEL% NEQ 0 GOTO end
      call %TestProgram% -C I_ADJUST_ADJUST_TI_CONTROL TI 10p
      IF %ERRORLEVEL% NEQ 0 GOTO end

210      REM --- MONITORED VALUES --- compares vent measurements to PTS2000 measure
      REM Vent measured PIP/PEEP Tolerance 8% + 2cmH2O

      call %TestProgram% -C I_COMPUTE_COMPUTED_VTI VTI 10p+10
215      IF %ERRORLEVEL% NEQ 0 GOTO end
      call %TestProgram% -C E_COMPUTE_COMPUTED_VTE VTI 10p+10

```

```

IF %ERRORLEVEL% NEQ 0 GOTO end

call %TestProgram% -C I_COMPUTE_MEASURE_R RATE 1
220 IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -C PMAX I_COMPUTE_COMPUTED_PIP/10 4p+2
IF %ERRORLEVEL% NEQ 0 GOTO end

call %TestProgram% -C PEEP E_COMPUTE_COMPUTED_PEEP/10 %peep_tolerance%
225 IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -C I_COMPUTE_MEASURE_TI/E_COMPUTE_MEASURE_TE IE 10p
IF %ERRORLEVEL% NEQ 0 GOTO end

230 REM ---- Alarm Checks: Check High Pressure Alarm is not set
call %TestProgram% -C HIGH_PRESSURE_ALARM I_ALARMS_ALARM_HIGH_PRESSURE
IF %ERRORLEVEL% NEQ 0 GOTO end

235 REM ---- Alarm Checks: Check Low Pressure Disconnection is not set
call %TestProgram% -C DISCONNECTION_ALARM I_ALARMS_ALARM_DISCONNECTION
IF %ERRORLEVEL% NEQ 0 GOTO end

240 REM ---- Starts ventilation ---
call %TestProgram% -R KEYBOARD_KEYBOARD_START_VENTIL_UP_EVENT 1
IF %ERRORLEVEL% NEQ 0 GOTO end

245 REM --- Perform 1 breath test after 20 breaths complete ---
call %TestProgram% -E %title% 20 1
IF %ERRORLEVEL% NEQ 0 GOTO end

250 REM --- VENTILATION STOP ---
call %TestProgram% -R KEYBOARD_KEYBOARD_STOP_VENTIL_EVENT 1
IF %ERRORLEVEL% NEQ 0 GOTO end

255 REM ----- Pressure AC Ped 1 Test -----

SET title="Mode=PRES, Pres=10cmH2O, Rate=15, Insp time=1.5s, Peep=0, Insp

260 REM modes: VOL=0 PRES=1 PSVT=2 VSIMV=3 CPAP=4
call %TestProgram% -R CONFIG_ADJUST_MODE 1
IF %ERRORLEVEL% NEQ 0 GOTO end

265 REM ----- Set Alarm Limits to OFF -----
call %TestProgram% -R ADJUST_LOW_VTI_NO_SELECT 1
IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -R ADJUST_HIGH_VTI_NO_SELECT 1
IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -R ADJUST_LOW_VTE_NO_SELECT 1
IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -R ADJUST_HIGH_VTE_NO_SELECT 1
IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -R CONFIG_FIO2_HIGH_SELECT 0
270 IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -R CONFIG_FIO2_LOW_SELECT 0
IF %ERRORLEVEL% NEQ 0 GOTO end

275 REM Turn OFF Barometric Compensation (0 = OFF)
call %TestProgram% -R CONFIG_ADJUST_BAROMETRIC_CORRECTOR 0
IF %ERRORLEVEL% NEQ 0 GOTO end

280 REM set PEEP to 0
call %TestProgram% -R ADJUST_PEEP_NO_SELECT 0
IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -R ADJUST_ADJUST_PEEP 0
IF %ERRORLEVEL% NEQ 0 GOTO end

285 REM RAMP 0
call %TestProgram% -R ADJUST_ADJUST_RAMP 0
IF %ERRORLEVEL% NEQ 0 GOTO end

REM Insp Time

```

```

290      call %TestProgram% -R ADJUST_ADJUST_TI_CONTROL 1500

      REM allow insp sens trigger setting to 5
      call %TestProgram% -R ADJUST_ADJUST_INSP_SENS 6
      IF %ERRORLEVEL% NEQ 0 GOTO end

295      REM Set R-Rate
      call %TestProgram% -R ADJUST_ADJUST_CONTROL_R 15
      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM Set Pressure Control (X 10)
300      call %TestProgram% -R ADJUST_ADJUST_P_CONTROL 100
      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM ----- REFRESH VENTILATOR SCREEN -----
305      REM (cursor needs to be on the first setting at the top of the VOL menu)
      call %TestProgram% -R KEYBOARD_KEYBOARD_UP_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
310      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
315      call %TestProgram% -R KEYBOARD_KEYBOARD_DOWN_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM -----

      REM --- SETTINGS COMPARISONS--- compares vent settings to PTS2000 measurem

320      SET peep_tolerance=1.5
      call %TestProgram% -C I_ADJUST_ADJUST_P_CONTROL+I_ADJUST_ADJUST_PEEP/1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      call %TestProgram% -C I_ADJUST_ADJUST_CONTROL_R RATE 1
325      IF %ERRORLEVEL% NEQ 0 GOTO end
      call %TestProgram% -C I_ADJUST_ADJUST_PEEP/10 PEEP %peep_tolerance%
      IF %ERRORLEVEL% NEQ 0 GOTO end
      call %TestProgram% -C I_ADJUST_ADJUST_TI_CONTROL TI 10p
      IF %ERRORLEVEL% NEQ 0 GOTO end

330

      REM --- MONITORED VALUES --- compares Vent measurements to PTS2000 measure
      REM Vent measured PIP/PEEP Tolerance 8% + 2cmH2O

335      call %TestProgram% -C I_COMPUTE_MEASURE_R RATE 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      call %TestProgram% -C PMAX I_COMPUTE_COMPUTED_PIP/10 4p+2
      IF %ERRORLEVEL% NEQ 0 GOTO end

340      call %TestProgram% -C PEEP E_COMPUTE_COMPUTED_PEEP/10 %peep_tolerance%
      IF %ERRORLEVEL% NEQ 0 GOTO end
      call %TestProgram% -C I_COMPUTE_MEASURE_TI/E_COMPUTE_MEASURE_TE IE 10p
      IF %ERRORLEVEL% NEQ 0 GOTO end

345

      REM ---- Alarm Checks: Check High Pressure Alarm is not set
      call %TestProgram% -C HIGH_PRESSURE_ALARM I_ALARMS_ALARM_HIGH_PRESSURE
      IF %ERRORLEVEL% NEQ 0 GOTO end

350      REM ---- Alarm Checks: Check Low Pressure Disconnection is not set
      call %TestProgram% -C DISCONNECTION_ALARM I_ALARMS_ALARM_DISCONNECTION
      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM ---- Starts ventilation ---
355      call %TestProgram% -R KEYBOARD_KEYBOARD_START_VENTIL_UP_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM --- Perform 1 breath test after 20 breaths complete ---
360      call %TestProgram% -E %title% 20 1
      IF %ERRORLEVEL% NEQ 0 GOTO end

```

```

365      REM ----- VENTILATION STOP -----
      call %TestProgram% -R KEYBOARD_KEYBOARD_STOP_VENTIL_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end

370      REM ----- Pressure AC Ped 2 Test -----

      SET title="Mode=PRES, Pres=5cmH2O, Rate=15, Insp time=1.5s, Peep=5, Insp S

      REM modes: VOL=0 PRES=1 PSVT=2 VSIMV=3 CPAP=4
      call %TestProgram% -R CONFIG_ADJUST_MODE 1
      IF %ERRORLEVEL% NEQ 0 GOTO end

375      REM Turn OFF Barometric Compensation (0 = OFF)
      call %TestProgram% -R CONFIG_ADJUST_BAROMETRIC_CORRECTOR 0
      IF %ERRORLEVEL% NEQ 0 GOTO end

380      REM set PEEP to adj mode
      call %TestProgram% -R ADJUST_PEEP_NO_SELECT 0
      IF %ERRORLEVEL% NEQ 0 GOTO end
      call %TestProgram% -R ADJUST_ADJUST_PEEP 50
      IF %ERRORLEVEL% NEQ 0 GOTO end

385      REM RAMP 0
      call %TestProgram% -R ADJUST_ADJUST_RAMP 0
      IF %ERRORLEVEL% NEQ 0 GOTO end

390      REM Insp Time
      call %TestProgram% -R ADJUST_ADJUST_TI_CONTROL 1500

      REM allow insp sens trigger setting to 5
      call %TestProgram% -R ADJUST_ADJUST_INSP_SENS 6
395      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM Set R-Rate
      call %TestProgram% -R ADJUST_ADJUST_CONTROL_R 15
      IF %ERRORLEVEL% NEQ 0 GOTO end

400      REM Set Pressure Control (X 10)
      call %TestProgram% -R ADJUST_ADJUST_P_CONTROL 50
      IF %ERRORLEVEL% NEQ 0 GOTO end

405      REM ----- REFRESH VENTILATOR SCREEN -----
      REM (cursor needs to be on the first setting at the top of the main menu i

      call %TestProgram% -R KEYBOARD_KEYBOARD_UP_EVENT 1
410      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
415      call %TestProgram% -R KEYBOARD_KEYBOARD_DOWN_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
420      REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
425      call %TestProgram% -R KEYBOARD_KEYBOARD_UP_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM
430      call %TestProgram% -R KEYBOARD_KEYBOARD_DOWN_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      REM -----

```

```

435 REM --- SETTINGS COMPARISONS--- compares vent settings to PTS2000 measurem
call %TestProgram% -C I_ADJUST_ADJUST_P_CONTROL+I_ADJUST_ADJUST_PEEP/1
IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -C I_ADJUST_ADJUST_CONTROL_R RATE 1
IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -C I_ADJUST_ADJUST_PEEP/10 PEEP 10p+1
440 IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -C I_ADJUST_ADJUST_TI_CONTROL TI 10p
IF %ERRORLEVEL% NEQ 0 GOTO end

445 REM --- MONITORED VALUES --- compares Vent measurements to PTS2000 measure
REM Vent measured PIP/PEEP Tolerance 8% + 2cmH2O

call %TestProgram% -C I_COMPUTE_MEASURE_R RATE 1
IF %ERRORLEVEL% NEQ 0 GOTO end
450 call %TestProgram% -C PMAX I_COMPUTE_COMPUTED_PIP/10 4p+2
IF %ERRORLEVEL% NEQ 0 GOTO end

call %TestProgram% -C PEEP E_COMPUTE_COMPUTED_PEEP/10 4p+2
IF %ERRORLEVEL% NEQ 0 GOTO end
455 call %TestProgram% -C I_COMPUTE_MEASURE_TI/E_COMPUTE_MEASURE_TE IE 10p
IF %ERRORLEVEL% NEQ 0 GOTO end

460 REM ---- Alarm Checks: Check High Pressure Alarm is not set
call %TestProgram% -C HIGH_PRESSURE_ALARM I_ALARMS_ALARM_HIGH_PRESSURE
IF %ERRORLEVEL% NEQ 0 GOTO end

REM ---- Alarm Checks: Check Low Pressure Disconnection is not set
call %TestProgram% -C DISCONNECTION_ALARM I_ALARMS_ALARM_DISCONNECTION
465 IF %ERRORLEVEL% NEQ 0 GOTO end

REM ---- Starts ventilation ---
call %TestProgram% -R KEYBOARD_KEYBOARD_START_VENTIL_UP_EVENT 1
470 IF %ERRORLEVEL% NEQ 0 GOTO end

REM --- Perform 1 breath test after 20 breaths complete ---
call %TestProgram% -E %title% 20 1
IF %ERRORLEVEL% NEQ 0 GOTO end
475

REM ----- VENTILATION STOP -----
call %TestProgram% -R KEYBOARD_KEYBOARD_STOP_VENTIL_EVENT 1
IF %ERRORLEVEL% NEQ 0 GOTO end

480

REM ----- Safety Pressure Test -----

485 SET title="Mode= Safety Pressure Test, Vol=400mL, Rate=15, Insp Time=1.5s,

REM modes: VOL=0 PRES=1 PSVT=2 VSIMV=3 CPAP=4
call %TestProgram% -R CONFIG_ADJUST_MODE 0
IF %ERRORLEVEL% NEQ 0 GOTO end
490

REM Set PIP alarm limits 20-900 (note: 20 equals 2) Set High limit = 90
call %TestProgram% -R ADJUST_ADJUST_LOW_PIP 20
IF %ERRORLEVEL% NEQ 0 GOTO end
495 call %TestProgram% -R ADJUST_ADJUST_HIGH_PIP 900
IF %ERRORLEVEL% NEQ 0 GOTO end

REM Turn OFF Barometric Compensation (0 = OFF)
call %TestProgram% -R CONFIG_ADJUST_BAROMETRIC_CORRECTOR 0
500 IF %ERRORLEVEL% NEQ 0 GOTO end

REM set PEEP to 0
call %TestProgram% -R ADJUST_PEEP_NO_SELECT 0
IF %ERRORLEVEL% NEQ 0 GOTO end

```

```

505      call %TestProgram% -R ADJUST_ADJUST_PEEP 0
      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM RAMP 0 = Slope, RAMP 1 = Rectangle
      call %TestProgram% -R ADJUST_ADJUST_RAMP 0
510      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM Insp Time
      call %TestProgram% -R ADJUST_ADJUST_TI_CONTROL 1500

515      REM allow insp sens trigger setting to 5
      call %TestProgram% -R ADJUST_ADJUST_INSP_SENS 6
      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM Set R-Rate
520      call %TestProgram% -R ADJUST_ADJUST_CONTROL_R 15
      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM Set volume (mL)
      call %TestProgram% -R ADJUST_ADJUST_VOL_CONTROL 400
525      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM --- REFRESH VENTILATOR SCREEN --- (cursor needs to be on the first set

530      call %TestProgram% -R KEYBOARD_KEYBOARD_UP_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      rem FOR /L %p IN (0,1,1000) DO REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
535      rem FOR /L %p IN (0,1,1000) DO REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
      rem FOR /L %p IN (0,1,1000) DO REM
      call %TestProgram% -R KEYBOARD_KEYBOARD_DOWN_EVENT 1
540      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM --- Safety Pressure Limit Test: Check Pressure to be 70-94cmH2O
      REM (TestPressure=82 +/- 12 [70-94cmH2O])
      SET TestPressure=82
545      call %TestProgram% -C %TestPressure% PMAX 12
      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM ----- Start ventilation -----
      call %TestProgram% -R KEYBOARD_KEYBOARD_START_VENTIL_UP_EVENT 1
550      IF %ERRORLEVEL% NEQ 0 GOTO end

      REM --- Perform 1 breath test after 10 breath complete ---
      call %TestProgram% -E %title% 5 1
      IF %ERRORLEVEL% NEQ 0 GOTO end
555

      REM ----- VENTILATION STOP -----
      call %TestProgram% -R KEYBOARD_KEYBOARD_STOP_VENTIL_EVENT 1
      IF %ERRORLEVEL% NEQ 0 GOTO end

560

      REM ----- High Pressure Alarm Test -----

565      SET title="Mode= High Pressure Alarm Test, Vol=400mL, Rate=15, Insp Time=1

      REM modes: VOL=0 PRES=1 PSVT=2 VSIMV=3 CPAP=4
      call %TestProgram% -R CONFIG_ADJUST_MODE 0
      IF %ERRORLEVEL% NEQ 0 GOTO end
570

      REM Set PIP alarm limits 20-600 (note: 20 equals 2) High limit = 60
      call %TestProgram% -R ADJUST_ADJUST_LOW_PIP 20
      IF %ERRORLEVEL% NEQ 0 GOTO end
575      call %TestProgram% -R ADJUST_ADJUST_HIGH_PIP 600
      IF %ERRORLEVEL% NEQ 0 GOTO end

```



```

    REM Turn OFF Barometric Compensation (0 = OFF)
    call %TestProgram% -R CONFIG_ADJUST_BAROMETRIC_CORRECTOR 0
580 IF %ERRORLEVEL% NEQ 0 GOTO end

    REM set PEEP to 0
    call %TestProgram% -R ADJUST_PEEP_NO_SELECT 0
    IF %ERRORLEVEL% NEQ 0 GOTO end
585 call %TestProgram% -R ADJUST_ADJUST_PEEP 0
    IF %ERRORLEVEL% NEQ 0 GOTO end

    REM RAMP 0 = Slope, RAMP 1 = Rectangle
    call %TestProgram% -R ADJUST_ADJUST_RAMP 0
590 IF %ERRORLEVEL% NEQ 0 GOTO end

    REM Insp Time
    call %TestProgram% -R ADJUST_ADJUST_TI_CONTROL 1500

595 REM allow insp sens trigger setting to 5
    call %TestProgram% -R ADJUST_ADJUST_INSP_SENS 6
    IF %ERRORLEVEL% NEQ 0 GOTO end

    REM Set R-Rate
    call %TestProgram% -R ADJUST_ADJUST_CONTROL_R 15
600 IF %ERRORLEVEL% NEQ 0 GOTO end

    REM Set volume (mL)
    call %TestProgram% -R ADJUST_ADJUST_VOL_CONTROL 400
605 IF %ERRORLEVEL% NEQ 0 GOTO end

    REM --- REFRESH VENTILATOR SCREEN --- (cursor needs to be on the first set

610 call %TestProgram% -R KEYBOARD_KEYBOARD_UP_EVENT 1
    IF %ERRORLEVEL% NEQ 0 GOTO end
    rem FOR /L %p IN (0,1,1000) DO REM
    call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
    IF %ERRORLEVEL% NEQ 0 GOTO end
615 rem FOR /L %p IN (0,1,1000) DO REM
    call %TestProgram% -R KEYBOARD_KEYBOARD_VALID_EVENT 1
    IF %ERRORLEVEL% NEQ 0 GOTO end
    rem FOR /L %p IN (0,1,1000) DO REM
    call %TestProgram% -R KEYBOARD_KEYBOARD_DOWN_EVENT 1
620 IF %ERRORLEVEL% NEQ 0 GOTO end

    REM ---- High Pressure Alarm Triggering: Check for High Pressure Alarm Eve
    call %TestProgram% -C 3 I_ALARMS_ALARM_HIGH_PRESSURE 0
625 IF %ERRORLEVEL% NEQ 0 GOTO end

    REM ---- Alarm Checks: Check Low Pressure Disconnection is not set
    call %TestProgram% -C DISCONNECTION_ALARM I_ALARMS_ALARM_DISCONNECTION
    IF %ERRORLEVEL% NEQ 0 GOTO end
630

    REM ---- Alarm Checks: Check Proximal Pressure Disconnection is not set
    call %TestProgram% -C 0 I_ALARMS_ALARM_PROXI_DISCONNECTION 1
    IF %ERRORLEVEL% NEQ 0 GOTO end

635 REM ----- Start ventilation -----
    call %TestProgram% -R KEYBOARD_KEYBOARD_START_VENTIL_UP_EVENT 1
    IF %ERRORLEVEL% NEQ 0 GOTO end

    REM --- Perform 1 breath test after 5 breath complete ---
640 call %TestProgram% -E %title% 5 1
    IF %ERRORLEVEL% NEQ 0 GOTO end

    REM ----- VENTILATION STOP -----
    call %TestProgram% -R KEYBOARD_KEYBOARD_STOP_VENTIL_EVENT 1
645 IF %ERRORLEVEL% NEQ 0 GOTO end

```

REM ----- Set Default Vent Settings -----

```

650 REM ----- in Volume mode -----

REM modes: VOL=0 PRES=1 PSVT=2 VSIMV=3 CPAP=4
call %TestProgram% -R CONFIG_ADJUST_MODE 0
IF %ERRORLEVEL% NEQ 0 GOTO end

655 REM set low Vte alarm limit to 300
call %TestProgram% -R ADJUST_LOW_VTE_NO_SELECT 0
IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -R ADJUST_ADJUST_LOW_VTE 300
660 IF %ERRORLEVEL% NEQ 0 GOTO end

REM set high Vte alarm limit to 1000
call %TestProgram% -R ADJUST_HIGH_VTE_NO_SELECT 0
IF %ERRORLEVEL% NEQ 0 GOTO end
665 call %TestProgram% -R ADJUST_ADJUST_HIGH_VTE 1000
IF %ERRORLEVEL% NEQ 0 GOTO end

REM set PEEP to 0
call %TestProgram% -R ADJUST_PEEP_NO_SELECT 0
670 IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -R ADJUST_ADJUST_PEEP 0
IF %ERRORLEVEL% NEQ 0 GOTO end

REM RAMP 0 = Slope, RAMP 1 = Rectangle
675 call %TestProgram% -R ADJUST_ADJUST_RAMP 0
IF %ERRORLEVEL% NEQ 0 GOTO end

REM Insp Time
call %TestProgram% -R ADJUST_ADJUST_TI_CONTROL 1300
680

REM allow insp sens trigger setting to 5
call %TestProgram% -R ADJUST_ADJUST_INSP_SENS 3
IF %ERRORLEVEL% NEQ 0 GOTO end
685

REM Set R-Rate
call %TestProgram% -R ADJUST_ADJUST_CONTROL_R 13
IF %ERRORLEVEL% NEQ 0 GOTO end

690 REM Set volume (mL)
call %TestProgram% -R ADJUST_ADJUST_VOL_CONTROL 500
IF %ERRORLEVEL% NEQ 0 GOTO end

REM ----- set pres mode -----
695 REM Set default settings: Pres A/C, P Control = 15, Peep = OFF, Rate = 13
REM Set I/E = 1:2.0, rise time = 2

REM modes: VOL=0 PRES=1 PSVT=2 VSIMV=3 CPAP=4
call %TestProgram% -R CONFIG_ADJUST_MODE 1
700 IF %ERRORLEVEL% NEQ 0 GOTO end

REM Turn ON Barometric Compensation (1 = ON)
call %TestProgram% -R CONFIG_ADJUST_BAROMETRIC_CORRECTOR 1
IF %ERRORLEVEL% NEQ 0 GOTO end
705

REM set low Vte alarm limit to 300
call %TestProgram% -R ADJUST_LOW_VTE_NO_SELECT 0
IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -R ADJUST_ADJUST_LOW_VTE 300
710 IF %ERRORLEVEL% NEQ 0 GOTO end

REM set high Vte alarm limit to 1000
call %TestProgram% -R ADJUST_HIGH_VTE_NO_SELECT 0
IF %ERRORLEVEL% NEQ 0 GOTO end
715 call %TestProgram% -R ADJUST_ADJUST_HIGH_VTE 1000
IF %ERRORLEVEL% NEQ 0 GOTO end

REM set low Vti alarm limit to 300
call %TestProgram% -R ADJUST_LOW_VTI_NO_SELECT 0
720 IF %ERRORLEVEL% NEQ 0 GOTO end

```

```

call %TestProgram% -R ADJUST_ADJUST_LOW_VTI 300
IF %ERRORLEVEL% NEQ 0 GOTO end

REM set high Vti alarm limit to 2000
725 call %TestProgram% -R ADJUST_HIGH_VTI_NO_SELECT 0
IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -R ADJUST_ADJUST_HIGH_VTI 2000
IF %ERRORLEVEL% NEQ 0 GOTO end

730 REM set High R = OFF
call %TestProgram% -R ADJUST_HIGH_R_NO_SELECT 1
IF %ERRORLEVEL% NEQ 0 GOTO end

REM set PEEP to 0
735 call %TestProgram% -R ADJUST_PEEP_NO_SELECT 0
IF %ERRORLEVEL% NEQ 0 GOTO end
call %TestProgram% -R ADJUST_ADJUST_PEEP 0
IF %ERRORLEVEL% NEQ 0 GOTO end

740 REM Insp Time
call %TestProgram% -R ADJUST_ADJUST_TI_CONTROL 1300

REM allow insp sens trigger setting
745 call %TestProgram% -R ADJUST_ADJUST_INSP_SENS 3
IF %ERRORLEVEL% NEQ 0 GOTO end

REM Set R-Rate
750 call %TestProgram% -R ADJUST_ADJUST_CONTROL_R 13
IF %ERRORLEVEL% NEQ 0 GOTO end

REM Set Pressure Control
755 call %TestProgram% -R ADJUST_ADJUST_P_CONTROL 150
IF %ERRORLEVEL% NEQ 0 GOTO end

REM ----- Relative pressure 'Off' -----
call %TestProgram% -R CONFIG_PRESSURE_SUPPORT_RELATIVE 0
760 IF %ERRORLEVEL% NEQ 0 GOTO end

REM ----- Ped Circuit 'Off' -----
call %TestProgram% -R CONFIG_ADJUST_PEDIATRIC_CIRCUIT 0
IF %ERRORLEVEL% NEQ 0 GOTO end
765 goto final

:end

PAUSE

770 REM ----- ERROR DETAILS -----
IF %ERRORLEVEL% == 1 goto syntax_error

IF %ERRORLEVEL% == 2 goto pts_error

775 IF %ERRORLEVEL% == 3 goto ventil_error

IF %ERRORLEVEL% == 4 goto oxy_error

780 IF %ERRORLEVEL% == 5 goto internal_error

goto final

785 :syntax_error

echo *****
echo *
echo * Syntax error *
790 echo *
echo *****

```

```
        goto final_error
795      :pts_error

        echo *****
        echo *
        echo *          PTS com error          *
800      echo *
        echo *****

        goto final_error

805      :ventil_error

        echo *****
        echo *
        echo *          Ventil com error        *
810      echo *
        echo *****

        goto final_error

815      :oxy_error

        echo *****
        echo *
        echo *          Oxygraph com error      *
820      echo *
        echo *****

        goto final_error

825      :internal_error

        echo *****
        echo *
        echo *          Internal error            *
830      echo *
        echo *****

        goto final_error

835      :final_error

        pause
        exit ERRORLEVEL

840      :final

        REM ----- Short Report generation -----

        call %GenerateProgram% -simple
845      exit 0
```