

REVISIONS	
REV	DESCRIPTION
F	RELEASE/CHANGE PER ECO-R179973

Test Protocol

CPU Board PB540, PB560

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 1 of 100	REV F
---	---	--	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

This document may not be copied, disclosed, or used in whole or in part without the consent of COVIDIEN.

© 2008 COVIDIEN All Rights Reserved.

Table of Contents

1.0	INTRODUCTION	4
1.1	PURPOSE	4
1.2	REVISION HISTORY	4
1.3	SCOPE	4
1.4	REFERENCE DOCUMENTS	4
1.5	ROLES AND RESPONSIBILITIES	4
1.6	ISSUE TRACKING	5
1.7	TEST EQUIPMENT	5
1.8	TEST UNITS	6
2.0	ACCEPTANCE CRITERIA	6
3.0	SAMPLE SIZE RATIONALE	7
3.1	VISUAL	7
3.2	POTENTIAL VARIATION	7
3.3	INSIGNIFICANT VARIATION	7
4.0	TEST COVERAGE	10
5.0	TEST PROCEDURES	11
5.1	24 VOLT SUPPLY MONITOR	11
5.2	POWER FAILURE INDICATOR	12
5.3	INRUSH CURRENT	13
5.4	CPU BOARD POWER CONSUMPTION	14
5.5	3.3 VOLT SUPPLY	15
5.6	5 VOLT SUPPLY	16
5.7	10 VOLT REFERENCE	17
5.8	5 VOLT REFERENCE	18
5.9	3.3 VOLT MONITOR	19
5.10	5 VOLT REFERENCE MONITOR	20
5.11	10 VOLT REFERENCE MONITOR	21
5.12	EXHALATION VALVE INTERFACE – PWM	22
5.13	EXHALATION VALVE INTERFACE – CURRENT	23
5.14	TURBINE INTERFACE – SPEED CONTROL PWM	24
5.15	TURBINE INTERFACE – BRAKE PWM	25
5.16	TURBINE INTERFACE – SPEED MEASUREMENT	26
5.17	TURBINE INTERFACE – ENABLE SIGNAL	27
5.18	TURBINE INTERFACE – TEMPERATURE MEASUREMENT	28
5.19	POWER SUPPLY COMMUNICATION INTERFACE – COMMUNICATION SPEED μ C	30
5.20	POWER SUPPLY COMMUNICATION INTERFACE – USB COMMUNICATION SPEED (PB560, PB520, FUTURE OPTION)	31
5.21	POWER SUPPLY COMMUNICATION INTERFACE – USB BUSY SIGNAL (PB560, PB520, FUTURE OPTION)	32
5.22	POWER SUPPLY COMMUNICATION INTERFACE – USB RESET SIGNAL	33
5.23	POWER SUPPLY COMMUNICATION INTERFACE – POWER INDICATORS	34
5.24	POWER SUPPLY COMMUNICATION INTERFACE – DATA TRANSFER SPEED	36
5.25	POWER SUPPLY COMMUNICATION INTERFACE – SOFTWARE DOWNLOAD	37
5.26	POWER SUPPLY COMMUNICATION INTERFACE – BOOTSTRAP MODE	38
5.27	POWER SUPPLY COMMUNICATION INTERFACE – BATTERY VOLTAGE MONITOR	39
5.28	BUZZER BOARD INTERFACE – PWM	40
5.29	BUZZER BOARD INTERFACE – LOGIC LEVELS	42
5.30	BUZZER BOARD INTERFACE – INHIBIT KEY	43
5.31	BUZZER BOARD INTERFACE – VOLTAGE MONITOR	44
5.32	REMOTE ALARM INTERFACE – N.C. RELAY CONTACTS	45
5.33	REMOTE ALARM INTERFACE – N.O. RELAY CONTACTS	46
5.34	KEYBOARD INTERFACE – “Down” KEY	47
5.35	KEYBOARD INTERFACE – “Up” KEY	48

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 2 of 100	REV F
---	---	--	-----------------

5.36	KEYBOARD INTERFACE – “VENTIL” KEY	49
5.37	KEYBOARD INTERFACE – “NAVIG” KEY	50
5.38	KEYBOARD INTERFACE – “VALID” KEY	51
5.39	KEYBOARD INTERFACE – “INHIB” KEY	52
5.40	KEYBOARD INTERFACE – BATTERY LED	53
5.41	KEYBOARD INTERFACE – DC PRESENCE LED	54
5.42	KEYBOARD INTERFACE – AC PRESENCE LED	55
5.43	KEYBOARD INTERFACE – OTHER LEDS	56
5.44	DISPLAY INTERFACE	57
5.45	O ₂ FLOW SENSOR INTERFACE (FUTURE OPTION)	58
5.46	O ₂ FLOW SENSOR INTERFACE (FUTURE OPTION)	59
5.47	O ₂ VALVE INTERFACE – PWM AND CURRENT DRIVE	60
5.48	SECURITY BACK-UP BUZZER INTERFACE – SECONDARY ALARM	62
5.49	SECURITY BACK-UP BUZZER INTERFACE – PRIMARY ALARM	63
5.50	RAM MEMORY	64
5.51	EVENT MEMORY	65
5.52	EVENT MEMORY – ERASE/PROGRAM CYCLES	66
5.53	VENTILATION SETTINGS MEMORY	67
5.54	VENTILATION SETTINGS MEMORY – ERASE/PROGRAM CYCLES	68
5.55	MONITORING MEMORY	69
5.56	MONITORING MEMORY – ERASE/PROGRAM CYCLES	70
5.57	WATCHDOG – POWER SUPERVISOR	71
5.58	WATCHDOG – TIMEOUT INTERVAL	74
5.59	WATCHDOG – PFI TRAP FUNCTION	75
5.60	CLOCK – FUNCTIONAL TEST	76
5.61	CLOCK – BATTERY TIMEKEEPING MODE	77
5.62	DEVICE MODEL ID	78
5.63	INSPIRATORY FLOW MEASUREMENT	79
5.64	EXHALATION FLOW MEASUREMENT	81
5.65	PROXIMAL PRESSURE MEASUREMENT	83
5.66	INTERNAL PRESSURE MEASUREMENT	85
5.67	EXPIRATORY VALVE PRESSURE MEASUREMENT	87
5.68	ATMOSPHERIC PRESSURE MEASUREMENT	89
5.69	O ₂ PRESSURE MEASUREMENT (FUTURE OPTION)	91
5.70	PCB – SIZE	93
5.71	PCB – MOUNTING	95
5.72	PCB – LABELING	96
5.73	PRINTED CIRCUIT BOARD FABRICATION STANDARDS	97
5.74	FIO ₂ MEASUREMENT (PB560, FUTURE OPTION)	98

TEST SUMMARY	99
---------------------------	-----------

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 3 of 100	REV F
---	---	--	-----------------

1.0 INTRODUCTION

1.1 Purpose

The purpose of this protocol is to test the CPU board used for PB540, PB560, PB520 and future ventilators.

1.2 Revision History

Revision	Author	Change Description
A		Incorporated review action items.
B		Clarification and/or correction of test steps.
C		Incorporate action Items from Design Verification Review
D		Add FiO2 testing to cover all the test for the PB560
E		Add USB Communication test to cover all the test for the PB560.
F		Corrected table in 3.3.4. Removed refs to 640. Removed test 75 as its not executed in this procedure – pointer to validation.

1.3 Scope

This protocol shall test the operational ranges and mechanical properties of the CPU board. Packaging standards, and tests which require an environmental chamber are excluded from this test protocol.

1.4 Reference Documents

1.4.1 Internal Documents

Ref	Part Number	Rev	Document Title
1	10024980	B	CPU Board Requirements
2	10035480	A	PRD PB560
3	10038637	A	PB560 Trace Matrix

1.5 Roles and Responsibilities

This section will define the roles and responsibilities required to release and execute this procedure.

R&D Engineer

The R&D Engineer will be responsible for the following activities:

- Generation of the Test Protocol
- Review of the Protocol
- Release of the Protocol, prior to execution

R&D Engineer / Test Engineer

An R&D Engineer or Test Engineer will be responsible for the following activities:

- Execution of the protocol
- Recording of test results
- Documentation of issues encountered
- Signing for completion

Review Team

The review team will consist of an R&D Engineer, the individual executing the protocol, a Quality Engineer, and an independent Reviewer.

The team will be responsible for the following activities:

- Review of the results
- Approval and acceptance of the execution results

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 4 of 100	REV F
---	---	--	-----------------

1.6 Issue Tracking

For the purpose of this protocol, issues will be logged in the following areas:

Issues encountered as part of development will be captured in Tracker, or as part of the action items during the review of the document.

Issues encountered during the execution of this protocol will be logged in the comments section of this document. Additional documents can be attached if required.

Issues discussed/encountered after execution will be captured as part of the review minutes or in Tracker.

1.7 Test Equipment

Test Equipment Description	Calibration ID	Calibration Due Date
PC with "ST10Flasher_BE" and "PICFlasher" software		
Digital power supply 0-36 VDC @ 4A		
Two multimeters with at least 4 digits of accuracy, and the ability to measure DC current.	Meter 1:	
	Meter 2:	
Digital oscilloscope.		
Current probe for oscilloscope.		
3020, 2588, 2228, 1924, and 1668 Ohm resistors, 1%, 0.1W	N/A	N/A
250 k Ohm variable resistor	N/A	N/A
Function generator able to generate a pulse with duty cycle adjustable from 50 to 80%, with frequency adjustable from 2 Hz to 1000 Hz, 0.1 % frequency accuracy.		
Flow meter from 0 to 1000 sccm with a \pm (3% of reading + 25 sccm) accuracy		
Pressure meter from -200 to +200 mbar \pm 0.1mbar		
Absolute pressure meter from 600 to 1100 mbar \pm 0.1mbar		
Sound Pressure Level meter		

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 5 of 100	REV F
---	---	--	-----------------

1.8 Test Units

Product Serial Number:	Unit 1:
	Unit 2:
	Unit 3:
	Unit 4:
	Unit 5:
	Unit 6:
Board Revision:	Unit 1:
	Unit 2:
	Unit 3:
	Unit 4:
	Unit 5:
	Unit 6:
Firmware version number:	Unit 1:
	Unit 2:
	Unit 3:
	Unit 4:
	Unit 5:
	Unit 6:
Software version number:	Unit 1:
	Unit 2:
	Unit 3:
	Unit 4:
	Unit 5:
	Unit 6:

2.0 ACCEPTANCE CRITERIA

The acceptance criteria will be included as part of the procedure to allow the individual executing the protocol to make the pass/fail determination.

For results that do not fall within the acceptance criteria, circle "Fail" and provide as much information in the comment section.

Some of the tests that have potential variation will have six (6) samples taken. The samples must satisfy a tolerance interval of 95% / 95% or have a rationale for not meeting the criteria.

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 6 of 100	REV F
---	---	--	-----------------

3.0 SAMPLE SIZE RATIONALE

The sample size will be based on the parameter type and verification method performed.

Initial testing will be performed to verify functionality by execution and analysis. Based on the test and parameter type, tests may be executed with multiple samples.

Testing will be conducted as specified in the criteria below

3.1 Visual

Parameters verified by visual examination, or analysis will be conducted with a sample size of 1.

3.1.1 Visual Exam

Verification of information contained on labels, confirmation of coloring, or verification of routing or hardware configuration is not expected to have significant variation

3.1.2 Datasheet / Specification Sheet Review

Verification done by review of data or specification sheets is not expected to have significant variation.

3.1.3 Design Analysis

Verification done through analysis is not expected to have significant variation.

3.2 Potential variation

A subset of the tests with parameters that have the potential for variation and have reproducible results will be tested with an initial sample of 6. Tests with variable input that have a variable range will not be used for data analysis.

Parameters that provide accuracy, peak values, or maximum thresholds that have variable (analog) data will be considered to have the potential of variation. The exception to this will be where a single, off the shelf, component influences the results (a parameter review may be sufficient to demonstrate adherence to the requirement).

3.2.1 Variable Measurements

- 3.2.1.1 Accuracy
- 3.2.1.2 Maximum level
- 3.2.1.3 Peak thresholds
- 3.2.1.4 Position

3.3 Insignificant variation

Testing of parameters, that are not expected to have significant variation (insignificant variation), will be conducted with three (3) samples as recommended in Sample Size Procedure (#10006201).

The following measurement types are not expected to provide significant variation.

3.3.1 Range

Verification that confirms operation or measurement of a range (non-accuracy) is not expected to have significant variation.

3.3.2 Single device/component verification

Verification of an output where a single off the shelf device (component) provides the resultant output is not expected to have significant variation. The specification of the device should demonstrate adherence to the requirement.

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 7 of 100	REV F
---	---	--	-----------------

3.3.3 Interconnect

Verification of an interconnect such as a trace or connector is not expected to have significant variation.

3.3.4 Digital

Measurements that provide binary, crystal controlled pulse widths, communication, Digital status, do not provide results where variation is expected or quantifiable.

3.3.4.1 Communication Interface

3.3.4.2 Digital I/O

3.3.4.3 Memory tests

3.3.4.4 Clocks

3.3.4.5 Pulse Width Modulation (PWM)

3.3.4.6 Initial Setting

3.3.4.7 Function/sequence operation

Test Procedure	Parameter/Test Type	Sample Size
24 Volt Supply Monitor	Range	3
Power Failure Indicator	Range	3
Inrush Current	Maximum limit	6
CPU Board Power Consumption	Peak threshold	6
3.3 Volt Supply	Single Device	3
5 Volt Supply	Single Device	3
10 Volt Reference	Single Device	3
5 Volt Reference	Single Device	3
3.3 Volt Monitor	Single Device	3
5 Volt Reference Monitor	Single Device	3
10 Volt Reference Monitor	Single Device	3
Exhalation Valve Interface – PWM	Digital	3
Exhalation Valve Interface – Current	Range	3
Turbine Interface – Speed Control PWM	Digital	3
Turbine Interface – Brake PWM	Digital	3
Turbine Interface – Speed Measurement	Accuracy	6
Turbine Interface – Enable Signal	Digital	3
Turbine Interface – Temperature Measurement	Range	3
Power Supply Communication Interface – Communication Speed uC	Digital	3
Power Supply Communication Interface – USB Communication Speed	Digital	3
Power Supply Communication Interface – USB Reset Signal	Digital	3
Power Supply Communication Interface – USB Busy Signal	Digital	3
Power Supply Communication Interface – Power Indicators	Digital	3
Power Supply Communication Interface – Data transfer speed	Digital	3
Power Supply Communication Interface – Device ID	Digital	3
Power Supply Communication Interface – Software Download	Digital	3
Power Supply Communication Interface – Bootstrap Mode	Digital	3

COVIDIEN
6135 Gunbarrel Avenue
Boulder, CO 80301

TITLE:
Test Protocol, CPU Board, PB540,
PB560

DOCUMENT NUMBER
10027876
SHEET 8 of 100

REV
F

Test Procedure	Parameter/Test Type	Sample Size
Power Supply Communication Interface – Battery Voltage Monitor	Range (accuracy)	6
Buzzer Board Interface – PWM	Digital	3
Buzzer Board Interface – Logic Levels	Digital	3
Buzzer Board Interface – Inhibit Key	Digital	3
Buzzer Board Interface – Voltage Monitor	Range (accuracy)	6
Remote Alarm Interface – N.C. Relay Contacts	Interconnect This test verifies resistance of the relay contacts similar to an interconnect	3
Remote Alarm Interface – N.O. Relay Contacts	Interconnect This verifies an open circuit (resistance of open relay).	3
Keyboard Interface – “Down” Key	Digital	3
Keyboard Interface – “Up” Key	Digital	3
Keyboard Interface – “Ventil” Key	Digital	3
Keyboard Interface – “Navig” Key	Digital	3
Keyboard Interface – “Valid” Key	Digital	3
Keyboard Interface – “Inhib” Key	Digital	3
Keyboard Interface – Battery LED	Digital	3
Keyboard Interface – DC Presence LED	Digital	3
Keyboard Interface – AC Presence LED	Digital	3
Keyboard Interface – Other LEDs	Digital	3
Display Interface	Digital	3
O2 Flow Sensor Interface (Future option)	Digital	3
O2 Flow Sensor Interface (Future option)	Digital	3
O2 Valve Interface – PWM and Current Drive	Range Digital - PWM	3
Security Back-up Buzzer Interface – Secondary Alarm	Digital	3
Security Back-up Buzzer Interface – Primary Alarm	Digital	3
RAM Memory	Digital	3
Event Memory	Digital	3
Event Memory – Erase/Program Cycles	Digital	3
Ventilation Settings Memory	Digital	3
Ventilation Settings Memory – Erase/Program Cycles	Digital	3
Monitoring Memory	Digital	3
Monitoring Memory – Erase/Program Cycles	Digital	3
Watchdog – Power Supervisor	Single device	3
Watchdog – Timeout interval	Digital	3
Watchdog – PFI Trap Function	Digital	3
Clock – Functional Test	Digital	3
Clock – Battery Timekeeping Mode	Digital	3
Device Model ID	Digital	3
Inspiratory Flow Measurement	Accuracy	6
Exhalation Flow Measurement	Accuracy	6
Proximal Pressure Measurement	Accuracy	6
Internal Pressure Measurement	Accuracy	6
Expiratory Valve Pressure Measurement	Accuracy	6
Atmospheric Pressure Measurement	Accuracy	6
O2 Pressure Measurement (Future option)	Accuracy	6
PCB – Size	Controlled on Drawing	3
PCB – Mounting	Controlled on Drawing	3
PCB – Labeling	Content	1

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 9 of 100	REV F
---	---	--	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

Test Procedure	Parameter/Test Type	Sample Size
Printed Circuit Board Fabrication Standards	Visual	1
FiO2 Measurement	Accuracy	6

4.0 TEST COVERAGE

The test cases in this procedure shall map to the component specification of the device.

For the initial release all test cases shall be executed.

A trace matrix will be included as part of the procedure, or attached, which shows 100% testing of the requirements contained in the component specification of the device, except packaging tests and system standards conformance which will be covered as part of system level testing.

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 10 of 100	REV F
---	---	---	-----------------

5.0 TEST PROCEDURES

5.1 24 Volt Supply Monitor

Test Number TCHWSCP B1; Referenced Requirement HWSCP B1.

5.1.1 Test Steps

- 1 - Disconnect the cable from the Power Management Board to the CPU board at the Power Management Board end.
- 2 - Connect the positive lead of a variable power supply, adjusted to 22 VDC, to J7-3 (+24VUTIL).
- 3 - Hold down the Inhibit key, then connect the negative lead of the power supply to J7-24 (Ground). Keep the key held down during startup to place the device into Maintenance Mode.
- 4 - Select "Measurements Check" from the menu.
- 5 - Verify that the displayed value of "24 V Check:" is 22 ± 0.5 VDC.
- 6 - Adjust the power supply to 24 VDC.
- 7 - Verify that the displayed value of "24 V Check:" is 24 ± 0.5 VDC.
- 8 - Adjust the power supply to 26 VDC.
- 9 - Verify that the displayed value of "24 V Check:" is 26 ± 0.5 VDC.

5.1.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
5 - Verify that the displayed "24 V Check:" voltage is 22 ± 0.5 VDC.	21.5 to 22.5	_____	Pass / Fail
7 - Verify that the displayed "24 V Check:" voltage is 24 ± 0.5 VDC.	23.5 to 24.5	_____	Pass / Fail
9 - Verify that the displayed "24 V Check:" voltage is 26 ± 0.5 VDC.	25.5 to 26.5	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

5.2 Power Failure Indicator

Test Number TCHWSCP2; Referenced Requirement HWSCP2.

5.2.1 Test Steps

- 1 - Disconnect the cable from the Power Management Board to the CPU board at the Power Management Board end.
- 2 - Connect the positive lead of a variable power supply, adjusted to 22 VDC, to J7-3 (24V UTIL).
- 3 - Hold down the Inhibit key, then connect the negative lead of the power supply to J7-24 (Ground). Keep the key held down during startup to place the device into Maintenance Mode.
- 4 - Select "Measurements Check" from the menu.
- 5 - Verify that the displayed voltage value of "Watchdog:" is 22 ± 0.5 VDC.
- 6 - Adjust the power supply to 24 VDC.
- 7 - Verify that the displayed voltage value of "Watchdog:" is 24 ± 0.5 VDC.
- 8 - Adjust the power supply to 26 VDC.
- 9 - Verify that the displayed voltage value of "Watchdog:" is 26 ± 0.5 VDC.

5.2.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
5 - Verify that the "Watchdog:" voltage value is 22 ± 0.5 VDC.	21.5 to 22.5	_____	Pass / Fail
7 - Verify that the "Watchdog:" voltage value is 24 ± 0.5 VDC.	23.5 to 24.5	_____	Pass / Fail
9 - Verify that the "Watchdog:" voltage value is 26 ± 0.5 VDC.	25.5 to 26.5	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 12 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

This document may not be copied, disclosed, or used in whole or in part without the consent of COVIDIEN.

© 2008 COVIDIEN All Rights Reserved.

5.3 Inrush Current

Test Number TCHWSCP3; Referenced Requirement HWSCP3.

5.3.1 Test Steps

- 1 - Disconnect the cable from the Power Management Board to the CPU board at the Power Management Board end.
- 2 - Use the regulated, Variable DC Power Supply (adjusted for 24 VDC) to provide power for the CPU Board. Set the power supply current limit to 3 A. Connect the positive power supply lead to J7-3 (24V UTIL), and the negative power supply lead to J7-24 (Ground).
- 3 - Use an oscilloscope with a current probe to monitor the current passing through the positive supply lead that is connected to J7-3.
- 4 - Turn on the power supply and capture the output of the current probe (the inrush current waveform) on the oscilloscope.
- 5 - Verify that the inrush current does not exceed 2 A.
- 6 - Verify that the inrush current is 100 ms or less in duration.

5.3.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
5 - Verify that the inrush current does not exceed 2 A.	< 2 A	_____	Pass / Fail
6 - Verify that the inrush current is 100 mS or less in duration.	< 100 mS	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 13 of 100	REV F
---	---	---	-----------------

5.4 CPU Board Power Consumption

Test Number TCHWSCP4; Referenced Requirement HWSCP4.

5.4.1 Test Steps

1 - Disconnect the cable from the Power Management Board to the CPU board at the Power Management Board end.

2 - Use the regulated, Variable DC Power Supply (adjusted for 24 VDC) to provide power for the CPU Board. Set the power supply current limit to 3 A. Connect the positive power supply lead to J7-3 (24V UTIL), and the negative power supply lead to J7-24 (Ground).

3 - Use an ammeter in-line with the positive supply lead (or the digital meter in the power supply, if available) to monitor the power supply current.

4 - Turn on the Power Supply and record the current to the CPU Board, in mA. (_____ mA)

5 - Calculate and record the power consumption (in milliwatts) as:

Power = 24 * (measured current in mA). (_____ mW)

6 - Verify that the power consumed from the 24 V Power Supply (calculated in the previous step) does not exceed 15000 mW (15 Watts.)

5.4.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
6 - Verify that the power consumed from the 24 V Power Supply does not exceed 15000 mW (15 Watts.)	< 15000 mW	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 14 of 100	REV F
---	---	---	-----------------

5.5 3.3 Volt Supply

Test Number TCHWSCP5; Referenced Requirement HWSCP8, HWSCP9 and HWSCP10.

5.5.1 Test Steps

- 1 - Referring to the specification sheet for the 3.3 Volt Regulator (IC26, LF33CDT) verify that the Output voltage is 3.3 ± 0.15 VDC at output currents between 0 and 400 mA DC, when provided an input voltage between 4.75 and 5.25 VDC.
- 2 - Connect a voltmeter to measure the DC voltage from the Output pin of IC26 to Ground.
- 3 - Verify that the voltmeter value is 3.3 ± 0.15 VDC.

5.5.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
1 - Verify that the regulator can provide an output voltage of 3.3 ± 0.15 VDC at output currents from 0 to 400 mA, when provided an input voltage from 4.75 to 5.25 VDC.	3.15 to 3.45 VDC	_____	Pass / Fail
3 - Verify that the voltmeter value is 3.3 ± 0.15 VDC.	3.15 to 3.45 VDC	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 15 of 100	REV F
---	---	---	-----------------

5.6 5 Volt Supply

Test Number TCHWSCP6; Referenced Requirement HWSCP5, HWSCP6 and HWSCP7.

5.6.1 Test Steps

- 1 - Referring to the specification sheet for the 5 Volt Regulator (IC29, TPS5430): Verify that the IC29 output voltage is 5 ± 0.25 VDC at output currents from 0 to 3 A, when provided an input voltage from 22.8 to 25.2 VDC.
- 2 - Measure the voltage across C106 with a voltmeter.
- 3 - Verify that the voltmeter value is 5 ± 0.25 VDC.

5.6.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
1 - Verify that the regulator can provide an output voltage of 5 ± 0.25 VDC at output currents from 0 to 3 A, when provided an input voltage from 22.8 to 25.2 VDC.	4.75 to 5.25 VDC	_____	Pass / Fail
3 - Verify that the voltmeter value is 5 ± 0.25 VDC.	4.75 to 5.25 VDC	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

5.7 10 Volt Reference

Test Number TCHWSCP B7; Referenced Requirement HWSCP B11, HWSCP B12 and HWSCP B13.

5.7.1 Test Steps

- 1 - Referring to the specification sheet for the 10 Volt Reference (IC23, LTS1236ACS8-10):
Verify that the IC23 output voltage is 10 ± 0.1 VDC at output currents from 0 to 10 mA, when provided an input voltage from 22.8 to 25.2 VDC.
- 2 - Measure the voltage across C85 with a voltmeter.
- 3 - Verify that the voltmeter value is 10 ± 0.1 VDC.

5.7.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
1 - Verify that the regulator can provide an output voltage of 10 ± 0.1 VDC at output currents from 0 to 10 mA, when provided an input voltage from 22.8 to 25.2 VDC.	9.9 to 10.1 VDC	_____	Pass / Fail
3 - Verify that the voltmeter value is 10 ± 0.1 VDC.	9.9 to 10.1 VDC	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 17 of 100	REV F
---	---	---	-----------------

5.8 5 Volt Reference

Test Number TCHWSCP8; Referenced Requirement HWSCP14, HWSCP15 and HWSCP16.

5.8.1 Test Steps

- 1 - Referring to the specification sheet for the 5 Volt Reference Regulator (IC21, LTS1236ACS8-5):
Verify that the IC21 output voltage is 5 ± 0.05 VDC at output currents from 0 to 10 mA, when provided an input voltage from 22.8 to 25.2 VDC.
- 2 - Measure the voltage across C81 with a voltmeter.
- 3 - Verify that the voltmeter value is 5 ± 0.05 VDC.

5.8.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
1 - Verify that the regulator can provide an output voltage of 5 ± 0.05 VDC at output currents from 0 to 10 mA, when provided an input voltage from 22.8 to 25.2 VDC.	4.95 to 5.05 VDC	_____	Pass / Fail
3 - Verify that the voltmeter value is 5 ± 0.05 VDC.	4.95 to 5.05 VDC	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

5.9 3.3 Volt Monitor

Test Number TCHWSCP B9; Referenced Requirement HWSCP B19 and HWSCP B20.

5.9.1 Test Steps

- 1 - Disconnect the cable between the CPU Board and the Power Management Board.
- 2 - Connect the positive lead of a programmable DC Power Supply (adjusted to 3.3 VDC) to IC36-3.
Connect the negative power supply lead to Ground (J7-24).
- 3 - Record the logic level value present on pin P2.5 (IC4-52, +3.3V-FAILURE).
- 4 - Verify that a HIGH logic level is present at P2.5 (IC4-52.).
- 5 - Observe P2.5 (IC4-52), while slowly decreasing the power supply voltage until pin P2.5 switches to a LOW logic level. Record the power supply voltage at which the signal on P2.5 went LOW.
- 6 - Verify that the power supply voltage is greater than or equal to 3 VDC.
- 7 - Adjust the programmable power supply to produce a repeating voltage which is at 3.3 VDC for 900 mS, and at 2.8 VDC for 100 mS of each cycle.
- 8 - Observe pin P2.5 (IC4-52) with an oscilloscope, and record the pulswidth of the LOW-going portion.
- 9 - Verify that the P2.5 remains at a LOW logic level for at least 140 mS of each cycle.

5.9.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that a HIGH logic level is present at P2.5 (IC4-52.)	HIGH logic level	_____	Pass / Fail
6 - Verify that the power supply voltage is greater than or equal to 3 VDC.	≥ 3 VDC	_____	Pass / Fail
9 - Verify that the P2.5 remains at a LOW logic level for at least 140 mS of each cycle.	LOW logic level	_____	Pass / Fail
9 - Verify that the P2.5 remains at a LOW logic level for at least 140 mS of each cycle.	≥ 140 mS	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 19 of 100	REV F
---	---	---	-----------------

5.10 5 Volt Reference Monitor

Test Number TCHWSCP10; Referenced Requirement HWSCP17 and HWSCP18.

5.10.1 Test Steps

- 1 - Disconnect the cable between the CPU Board and the Power Management Board.
- 2 - Connect the positive lead of a programmable DC Power Supply (set to 5 VDC) to IC24-3. Connect the negative power supply lead to Ground (J7-24). Monitor the power supply voltage with a voltmeter.
- 3 - Record the logic level value present on pin P2.6 (IC4-53, +5VREF-FAILURE).
- 4 - Verify that a HIGH logic level is present at P2.6 (IC4-53).
- 5 - Observe P2.6 (IC4-53), while slowly decreasing the power supply voltage until pin P2.6 switches to a LOW logic level. Record the power supply voltage at which the signal on P2.6 went LOW.
- 6 - Verify that the power supply voltage is greater than or equal to 4.5 VDC.
- 7 - Adjust the programmable power supply to produce a repeating voltage which is at 5 VDC for 900 mS, and at 4 VDC for 100 mS of each cycle.
- 8 - Observe pin P2.6 (IC4-53) with an oscilloscope, and record the pulsewidth of the LOW-going portion of the signal.
- 9 - Verify that the P2.6 remains at a LOW logic level for at least 140 mS of each cycle.

5.10.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that a HIGH logic level is present at P2.6 (IC4-53).	HIGH logic level	_____	Pass / Fail
6 - Verify that the power supply voltage is greater than or equal to 4.5 VDC.	≥ 4.5 VDC	_____	Pass / Fail
9 - Verify that P2.6 remains at a LOW logic level for at least 140 mS of each cycle.	LOW logic level	_____	Pass / Fail
9 - Verify that P2.6 remains at a LOW logic level for at least 140 mS of each cycle.	≥ 140 mS	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 20 of 100	REV F
---	---	---	-----------------

5.11 10 Volt Reference Monitor

Test Number TCHWSCP11; Referenced Requirement HWSCP21 and HWSCP22.

5.11.1 Test Steps

- 1 - Disconnect the cable between the CPU Board and the Power Management Board.
- 2 - Connect the positive lead of a 24 VDC power supply to J7-3, and the negative lead to J7-24 to provide power the CPU board.
- 3 - Unsolder the "top" end of R66 (the end furthest from the edge of the CPU Board) so that it does not make electrical contact with the CPU Board.
- 4 - Connect the positive lead of a programmable DC Power Supply (set to 10 VDC) to the "top" of R66. Connect the negative power supply lead to Ground (J7-24).
- 5 - Record the logic level value present on pin P2.7 (IC4-54, +10VREF-FAILURE).
- 6 - Verify that a HIGH logic level is present at P2.7 (IC4-54).
- 7 - Observe P2.7 (IC4-54), while slowly decreasing the power supply voltage until pin P2.7 switches to a LOW logic level. Record the power supply voltage at which the signal on P2.7 went LOW.
- 8 - Verify that the power supply voltage is greater than or equal to 8.9 VDC.
- 9 - Adjust the programmable power supply to produce a repeating voltage which is at 10 VDC for 900 mS, and at 8 VDC for 100 mS of each cycle.
- 10 - Observe pin P2.7 (IC4-54) with an oscilloscope, and record the pulsewidth of the LOW-going portion of the signal.
- 11 - Verify that the P2.7 remains at a LOW logic level for at least 140 mS of each cycle.
- 12 - Reattach R66 to the CPU board as it originally was.

5.11.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
6 - Verify that a HIGH logic level is present at P2.7 (IC4-54).	HIGH logic level	_____	Pass / Fail
8 - Verify that the power supply voltage is greater than or equal to 8.9 VDC.	≥ 8.9 VDC	_____	Pass / Fail
11 - Verify that P2.7 remains at a LOW logic level for at least 140 mS of each cycle.	LOW logic level	_____	Pass / Fail
11 - Verify that P2.7 remains at a LOW logic level for at least 140 mS of each cycle.	≥ 140 mS	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 21 of 100	REV F
---	---	---	-----------------

5.12 Exhalation Valve Interface – PWM

Test Number TCHWSCP12; Referenced Requirement HWSCP23, HWSCP25 and HWSCP26.

5.12.1 Test Steps

- 1 - Generate a 14.975 kHz PWM signal with a pulsewidth of 33.38 uS on microcontroller Port P7.1 (IC4-20) by setting the appropriate registers.
- 2 - Measure the frequency of the signal at G-T21 (CD-VALVE), using an oscilloscope.
- 3 - Verify that the frequency of the signal is 14.975 kHz \pm 10%.
- 4 - Generate a 14.975 kHz PWM signal with a pulsewidth of approximately 25 nS on pin P7.1 (IC4-20) by setting the appropriate registers.
- 5 - Measure the pulsewidth with the oscilloscope.
- 6 - Generate a 14.975 kHz PWM signal with a pulsewidth of approximately 50 nS on P7.1 by setting the appropriate registers.
- 7 - Measure the pulsewidth with the oscilloscope.
- 8 - Verify that the difference between first and second pulsewidth measurement is 25 \pm 5 nS.

5.12.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that the frequency of the signal is 14.975 kHz \pm 10%.	13.4775 to 16.4725	_____	Pass / Fail
5 - Measure the pulsewidth with the oscilloscope.	25 nS	_____	N/A
7 - Measure the pulsewidth with the oscilloscope.	50 nS	_____	N/A
8 - Verify that the difference between first and second pulsewidth measurement is 25 \pm 5 nS.	20 to 30 nS	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 22 of 100	REV F
---	---	---	-----------------

5.13 Exhalation Valve Interface – Current

Test Number TCHWSCP13; Referenced Requirement HWSCP24, HWSCP27 and HWSCP28.

5.13.1 Test Steps

- 1 - Unplug the J19 (Exhalation Valve) connector.
- 2 - Generate a HIGH logic level signal on P7.1 (IC4-20) by setting the appropriate microcontroller registers.
- 3 - Read the digital value of the MES-I-VALVE signal at P5.11 (IC4-40) by setting the appropriate microcontroller registers.
- 4 - Verify that the digital value is 0 ± 21 .
- 5 - Connect a 750 Ohm resistor, in series with an ammeter, between J19-1 and J19-2.
- 6 - Read the digital value of P5.11 (IC4-40) by setting the appropriate registers.
- 7 - Verify that the digital value is equal to $(\text{Ammeter value} * 20460) \pm 21$.
- 8 - Connect a 1500 Ohm resistor, in series with an ammeter, between J19-1 and J19-2.
- 9 - Read the digital value of P5.11 (IC4-40) by setting the appropriate registers.
- 10 - Verify that the digital value is equal to $(\text{Ammeter value} * 20460) \pm 21$.

5.13.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that the digital value is 0 ± 21 .	-21 to 21	_____	Pass / Fail
7 - Verify that the digital value is equal to $(\text{Ammeter value} * 20460) \pm 21$.	$(Av * 20460) \pm 21$	_____	Pass / Fail
10 - Verify that the digital value is equal to $(\text{Ammeter value} * 20460) \pm 21$.	$(Av * 20460) \pm 21$	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 23 of 100	REV F
---	---	---	-----------------

5.14 Turbine Interface – Speed Control PWM

Test Number TCHWSCP14; Referenced Requirement HWSCP29 and HWSCP30.

5.14.1 Test Steps

- 1 - Generate a 9.765 kHz PWM signal with a pulsewidth of 102.4 uS on P7.0 (IC4-19) by setting the appropriate registers.
- 2 - Observe the signal at J4-5 (SPEED-SETPOINT) with an oscilloscope:
- 3 - Verify that the frequency of the PWM signal is 9.765 kHz \pm 10%.
- 4 - Generate a 9.765 kHz PWM signal with a pulsewidth of approximately 50 nS on pin P7.0 (IC4-19) by setting the appropriate registers.
- 5 - Measure the pulsewidth with the oscilloscope.
- 6 - Generate a 9.765 kHz PWM signal with a pulsewidth of approximately 75 nS on pin P7.0 (IC4-19) by setting the appropriate registers.
- 7 - Measure the pulsewidth with the oscilloscope.
- 8 - Verify that the difference between first and second pulsewidth measurement is 25 \pm 5 nS.

5.14.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that the frequency of the PWM signal is 9.765 kHz \pm 10%.	8.7885 to 10.7415 kHz	_____	Pass / Fail
5 - Measure the pulsewidth with the oscilloscope.	50 nS	_____	N/A
7 - Measure the pulsewidth with the oscilloscope.	75 nS	_____	N/A
8 - Verify that the difference between first and second pulsewidth measurement is 25 \pm 5 nS.	20 to 30 nS	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 24 of 100	REV F
---	---	---	-----------------

5.15 Turbine Interface – Brake PWM

Test Number TCHWSCP15; Referenced Requirement HWSCP31 and HWSCP32.

5.15.1 Test Steps

- 1 - Generate a 15.060 kHz PWM signal with a pulsewidth of 33.2 uS on P2.1 (IC4-48) by setting the appropriate registers.
- 2 - Observe the signal at J4-4 (BRAKE/) with the oscilloscope, and measure its frequency.
- 3 - Verify that the frequency is 15.060 kHz \pm 10%.
- 4 - Generate a 15.060 kHz PWM signal with a pulsewidth of approximately 400 nS on P2.1 (IC4-48) by setting the appropriate registers.
- 5 - Measure pulsewidth with the oscilloscope.
- 6 - Generate a 15.060 kHz PWM signal with a pulsewidth of approximately 800 nS on P2.1 (IC4-48) by setting the appropriate registers.
- 7 - Measure pulsewidth with the oscilloscope.
- 8 - Verify that the difference between first and second pulsewidth measurement is 400 \pm 80 nS.

5.15.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that the frequency is 15.060 kHz \pm 10%.	13.554 to 16.566 kHz	_____	Pass / Fail
5 - Measure the pulsewidth with the oscilloscope.	400 nS	_____	N/A
7 - Measure the pulsewidth with the oscilloscope.	800 nS	_____	N/A
8 - Verify that the difference between first and second pulsewidth measurement is 400 \pm 80 nS.	320 to 480 nS	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 25 of 100	REV F
---	---	---	-----------------

5.16 Turbine Interface – Speed Measurement

Test Number TCHWSCP16; Referenced Requirement HWSCP33.

5.16.1 Test Steps

- 1 - Connect a function generator to J4-6 (SPEED-MEASURE) and apply a 0-to-5V 1 kHz square wave signal.
- 2 - Download approved “CPU Software” version to the CPU Board, and approved “Soft carte alim pour test CPU” software version to the Power Supply Management board.
- 3 - Start ventilator in “Setup” mode, then run the Maintenance page.
- 4 - Verify that the turbine speed measurement displayed on ventilator is 60,000 RPM \pm 1%
- 5 - Change function generator frequency to 500 Hz.
- 6 - Verify that the turbine speed measurement displayed on ventilator is 30,000 RPM \pm 1%
- 7 - Change function generator frequency to 16.67 Hz.
- 8 - Verify that the turbine speed measurement displayed on ventilator is 1,000 RPM \pm 1%

5.16.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that the turbine speed measurement displayed on ventilator is 60,000 RPM \pm 1%	59,400 to 60,600 RPM	_____	Pass / Fail
6 - Verify that the turbine speed measurement displayed on ventilator is 30,000 RPM \pm 1%	29,700 to 30,300 RPM	_____	Pass / Fail
8 - Verify that the turbine speed measurement displayed on ventilator is 1,000 RPM \pm 1%	990 to 1,010 RPM	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 26 of 100	REV F
---	---	---	-----------------

5.17 Turbine Interface – Enable Signal

Test Number TCHWSCP17; Referenced Requirement HWSCP34.

5.17.1 Test Steps

- 1 - Generate a LOW logic level at pin P2.3 (IC4-50, ENABLE-TURB) by setting the appropriate registers.
- 2 - Verify that a LOW logic level is present on J4-2 (ENABLE).
- 3 - Generate a HIGH logic level at pin P2.3 (IC4-50) by setting the appropriate registers.
- 4 - Verify that a HIGH logic level is present on J4-2.

5.17.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - Verify that a LOW logic level is present on J4-2 (ENABLE).	LOW logic level	_____	Pass / Fail
4 - Verify that a HIGH logic level is present on J4-2.	HIGH logic level	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 27 of 100	REV F
---	---	---	-----------------

5.18 Turbine Interface – Temperature Measurement

Test Number TCHWSCP18; Referenced Requirement HWSCP35.

5.18.1 Test Steps

- 1 - Connect a resistance decade box (or use suitable fixed resistors) between J4-10 (TURBINE-TEMP) and J4-9 (0VANA).
- 2 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Mgmt. Board.
- 3 - Start ventilator in "Partial Tests" mode, then run the Blower Interface test.
- 4 - Adjust resistance to 3020 Ohms.
- 5 - Verify that the blower temperature measurement displayed on the ventilator is 60 ± 1 °C.
- 6 - Adjust resistance to 2588 Ohms.
- 7 - Verify that the blower temperature displayed is 65 ± 1 °C.
- 8 - Adjust resistance to 2228 Ohms.
- 9 - Verify that the blower temperature displayed is 70 ± 1 °C.
- 10 - Adjust resistance to 1924 Ohms.
- 11 - Verify that the blower temperature displayed is 75 ± 1 °C.
- 12 - Adjust resistance to 1668 Ohms.
- 13 - Verify that the blower temperature displayed is 80 ± 1 °C.

5.18.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
5 - Verify that the blower temperature measurement displayed on the ventilator is 60 ± 1 °C.	59 to 61 °C	_____	Pass / Fail
7 - Verify that the blower temperature displayed is 65 ± 1 °C.	64 to 66 °C	_____	Pass / Fail
9 - Verify that the blower temperature displayed is 70 ± 1 °C.	69 to 71 °C	_____	Pass / Fail
11 - Verify that the blower temperature displayed is 75 ± 1 °C.	74 to 76 °C	_____	Pass / Fail
13 - Verify that the blower temperature displayed is 80 ± 1 °C.	79 to 81 °C	_____	Pass / Fail

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 28 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

Comments: _____

Test Operator Signature/Date: _____

5.19 Power Supply Communication Interface – Communication speed μ C

Test Number TCHWSCP B19; Referenced Requirement HWSCP B36.

5.19.1 Test Steps

- 1 - Observe pin J7-BR6 with an oscilloscope, and record the signal.
- 2 - Verify that the frequency of the signal is between 1.71Mbit/s and 1.89Mbit/s.
Warning: Measure only the frequency include into the burst.

5.19.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
(2) Communication shall have a transfer rate of 1.8Mbit/s +/-5%.	1.71Mbit/s to 1.89Mbit/s	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 30 of 100	REV F
---	---	---	-----------------

5.20 Power Supply Communication Interface – USB Communication speed (PB560, PB520, Future option)

Test Number TCHWSCP76; Referenced Requirement HWSCP117.

5.20.1 Test Steps

- 1 - Observe pin J7-BR5 with an oscilloscope, and record the signal.
- 2 - Verify that the frequency of the signal is between 1.71Mbit/s and 1.89Mbit/s.
Warning: Measure only the frequency include into the burst.

5.20.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
(2) Communication shall have a transfer rate of 1.8Mbit/s +/-5%.	1.71Mbit/s to 1.89Mbit/s	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 31 of 100	REV F
---	---	---	-----------------

5.21 Power Supply Communication Interface – USB Busy Signal (PB560, PB520, Future Option)

Test Number TCHWSCP B77; Referenced Requirement HWSCP B44.

5.21.1 Test Steps

- 1 - Generate a LOW logic level on J7-12 (USB-BUSY).
- 2 - Verify a LOW logic level at pin P2.10 (IC4-59) by setting the appropriate registers.
- 3 - Generate a HIGH logic level on J7-12 (USB-BUSY).
- 4 - Verify a HIGH logic level at pin P2.10 (IC4-59) by setting the appropriate registers.

5.21.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - Verify a LOW logic level.	LOW logic level	_____	Pass / Fail
4 - Verify a HIGH logic level.	HIGH logic level	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 32 of 100	REV F
---	---	---	-----------------

5.22 Power Supply Communication Interface – USB Reset Signal

Test Number TCHWSCP B20; Referenced Requirement HWSCP B43.

5.22.1 Test Steps

- 1 - Generate a LOW logic level at pin P3.5 (IC4-70) by setting the appropriate registers.
- 2 - Place a voltmeter on J7-11 (RESET-USB)
- 3 - Verify that the voltage is 0 ± 0.5 VDC.
- 4 - Generate a HIGH logic level at P3.5 by setting the appropriate registers.
- 5 -Verify that the voltage is 5 ± 0.5 VDC.

5.22.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that the voltage is 0 ± 0.5 VDC.	-0.5 to 0.5 VDC	_____	Pass / Fail
5 -Verify that the voltage is 5 ± 0.5 VDC.	4.5 to 5.5 VDC	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 33 of 100	REV F
---	---	---	-----------------

5.23 Power Supply Communication Interface – Power Indicators

Test Number TCHWSCP22; Referenced Requirement HWSCP37, HWSCP38 and HWSCP39.

5.23.1 Test Steps

- 1 - Use a regulated power supply to apply 0 VDC to J7-15 (LED-PRESENCE-AC).
- 2 - Verify that the AC Power LED is OFF.
- 3 - Read the logic level of pin P8.0 (IC4-9, AC-PRESENCE).
- 4 - Verify that a LOW logic level is present at P8.0 (IC4-9).
- 5 - Use a regulated power supply to apply 5 VDC to J7-15.
- 6 - Verify that the AC Power LED is ON.
- 7 - Read the logic level of pin P8.0 (IC4-9, AC-PRESENCE).
- 8 - Verify that a HIGH logic level is present at P8.0 (IC4-9).
- 9 - Use a regulated power supply to apply 0 VDC to J7-16 (LED-PRESENCE-DC).
- 10 - Verify that the DC Power LED is OFF.
- 11 - Use a regulated power supply to apply 5 VDC to J7-16.
- 12 - Verify that the DC Power LED is ON.
- 13 - Use a regulated power supply to apply 0 VDC to J7-17 (LED-BAT).
- 14 - Verify that the Internal Battery LED is OFF.
- 15 - Use a regulated power supply to apply 5 VDC to J7-17.
- 16 - Verify that the Internal Battery LED is ON.

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 34 of 100	REV F
---	---	---	-----------------

5.23.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - Verify that the AC Power LED is OFF.	LED off	_____	Pass / Fail
4 - Verify that a LOW logic level is present at P8.0 (IC4-9).	LOW logic level	_____	Pass / Fail
6 - Verify that the AC Power LED is ON.	LED on	_____	Pass / Fail
8 - Verify that a HIGH logic level is present at P8.0 (IC4-9).	HIGH logic level	_____	Pass / Fail
10 - Verify that the DC Power LED is OFF.	LED off	_____	Pass / Fail
12 - Verify that the DC Power LED is ON.	LED on	_____	Pass / Fail
14 - Verify that the Internal Battery LED is OFF.	LED off	_____	Pass / Fail
16 - Verify that the Internal Battery LED is ON.	LED on	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 35 of 100	REV F
---	---	---	-----------------

5.24 Power Supply Communication Interface – Data transfer speed

Test Number TCHWSCP B23; Referenced Requirement HWSCP B40.

5.24.1 Test Steps

- 1 – Connect the PB540 to a PC with the USB cable.
- 2 – Open “*ST10Flasher_BE*” software (this software has a communication rate set to 115,200bauds by default). Check “LEGENDAIR US” in the windows “S/N” on the right. Automatically the number “40971K” is written on the box bellow. Complete the S/N by the actual number of the PB540 tested – 6 numbers XXXXXX and click on the button “Ecriture S/N”. This function will write and read the S/N on the device, and will allow us to test the communication.
- 3 – Verify that on the bottom of the main windows of the software, the following phrase is written: “Numero de série 40971Kxxxxxx ecrit avec success”.
That's means that the software read the S/N in the device to ensure it has been written right.

5.24.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
(2) Communication shall be capable of data transfers with a baud rate of 115200.	Verify the phrase “Numero de série 40971Kxxxxxx ecrit avec success » on the PC software	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 36 of 100	REV F
---	---	---	-----------------

5.25 Power Supply Communication Interface – Software Download

Test Number TCHWSCP24; Referenced Requirement HWSCP41.

5.25.1 Test Steps

- 1 - Use "ST10Flasher" software to download an approved "CPU Test Software" version to the CPU Board at 256000 baud.
- 2 - Verify that the "Init Monitor", "Load File", "Erase Flash", and "Program Flash" fields have a status of "OK".

5.25.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - Verify that the "Init Monitor" field has a status of "OK".	"Init Monitor" OK	_____	Pass / Fail
2 - Verify that the "Load File" field has a status of "OK".	"Load File" OK	_____	Pass / Fail
2 - Verify that the "Erase Flash" field has a status of "OK".	"Erase Flash" OK	_____	Pass / Fail
2 - Verify that the "Program Flash" field has a status of "OK".	"Program Flash" OK	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 37 of 100	REV F
---	---	---	-----------------

5.26 Power Supply Communication Interface – Bootstrap Mode

Test Number TCHWSCP25; Referenced Requirement HWSCP42.

5.26.1 Test Steps

- 1 - Turn main power switch OFF.
- 2 - Use “ST10Flasher” software to enter bootstrap mode.
- 3 - Turn main power switch ON.
- 4 - Verify that the “Bootstrap” button changes to “Programming”, and the device screen stays cleared.

5.26.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that the bootstrap button changes to “Programming”.	“Programming”	_____	Pass / Fail
4 - Verify that the device screen stays cleared.	Screen clear	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 38 of 100	REV F
---	---	---	-----------------

5.27 Power Supply Communication Interface – Battery Voltage Monitor

Test Number TCHWSCP26; Referenced Requirement HWSCP45.

5.27.1 Test Steps

- 1 - Remove the battery from the device.
- 2 - Cover the positive battery terminal with tape, and reinstall the battery into the device.
- 3 - Hold down the "Inhibit" key while turning the device on to enter Maintenance mode.
- 4 - Select "Maintenance", then "Measurement Check", and then "Internal Battery Menu" from the menu screens.
- 5 - Use a variable power supply to apply 0 VDC to J7-14 (VBAT).
- 6 - Verify that the displayed Battery Voltage value is 0 ± 0.3 VDC.
- 7 - Adjust the variable power supply to apply 15 VDC to J7-14 (VBAT).
- 8 - Verify that the displayed Battery Voltage value is 15 ± 0.3 VDC.
- 9 - Adjust the variable power supply to apply 30 VDC to J7-14 (VBAT).
- 10 - Verify that the displayed Battery Voltage value is 30 ± 0.3 VDC.
- 11 - Remove the battery, remove tape from positive terminal, and reinstall battery.

5.27.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
6 - Verify that the displayed Battery Voltage value is 0 ± 0.3 VDC.	-0.3 to 0.3 VDC	_____	Pass / Fail
8 - Verify that the displayed Battery Voltage value is 15 ± 0.5 VDC.	14.7 to 15.3 VDC	_____	Pass / Fail
10 - Verify that the displayed Battery Voltage value is 30 ± 0.5 VDC.	29.7 to 30.3 VDC	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 39 of 100	REV F
---	---	---	-----------------

5.28 Buzzer Board Interface – PWM

Test Number TCHWSCP B27; Referenced Requirement HWSCP B46, HWSCP B47 and HWSCP B48.

5.28.1 Test Steps

- 1 - Observe the signal at J18-3 with an oscilloscope.
- 2 - Generate a 60 kHz PWM signal with a pulsewidth of 8.33 μ S on pin P7.3 (IC4-22, PWM-BUZ) by setting the appropriate registers.
- 3 - Verify with the oscilloscope that the modulated carrier frequency is 60 kHz \pm 10%.
- 4 - Generate an 880 Hz PWM signal with a pulsewidth of 0.57 mS on pin P7.4 (IC4-23) by setting the appropriate registers.
- 5 - Verify with the oscilloscope that the modulating signal frequency is 880 Hz \pm 10%.
- 6 - Generate a 1000 Hz PWM signal with a pulsewidth of 0.5 mS on pin P7.4 (IC4-23) by setting the appropriate registers.
- 7 - Verify with the oscilloscope that the modulating signal frequency is 1000 Hz \pm 10%.
- 8 - Generate a 60 kHz PWM signal with a pulsewidth of approximately 300 nS on pin P7.3 (IC4-22) by setting the appropriate registers.
- 9 - Measure pulsewidth with the oscilloscope.
- 10 - Generate a 60 kHz PWM signal with a pulsewidth of approximately 350 nS on pin P7.3 (IC4-22) by setting the appropriate registers.
- 11 - Measure pulsewidth with the oscilloscope.
- 12 - Verify that the difference between step 9 and step 11 pulsewidth measurements is 50 \pm 10ns.

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 40 of 100	REV F
---	---	---	-----------------

5.28.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify with the oscilloscope that the modulated carrier frequency is 60 kHz \pm 10%.	54 to 66 kHz	_____	Pass / Fail
5 - Verify with the oscilloscope that the modulating signal frequency is 880 Hz \pm 10%.	792 to 968 Hz	_____	Pass / Fail
7 - Verify with the oscilloscope that the modulating signal frequency is 1000 Hz \pm 10%.	900 to 1100 Hz	_____	Pass / Fail
9 - Measure pulsewidth with the oscilloscope.	N/A	_____	N/A
11 - Measure pulsewidth with the oscilloscope.	N/A	_____	N/A
12 - Verify that the difference between step 9 and step 11 pulsewidth measurements is 50 \pm 10ns.	40 to 60 ns	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 41 of 100	REV F
---	---	---	-----------------

5.29 Buzzer Board Interface – Logic Levels

Test Number TCHWSCP B28; Referenced Requirement HWSCP B49.

5.29.1 Test Steps

- 1 - Use an oscilloscope to capture and measure the voltage levels on J18-4 during startup.
- 2 - Verify that the LOW logic level voltage at J18-4 is 0 ± 0.5 VDC.
- 3 - Verify that the HIGH logic level voltage at J18-4 is 5 ± 0.5 VDC.

5.29.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - Verify that the LOW logic level voltage at J18-4 is 0 ± 0.5 VDC.	-0.5 to 0.5 VDC	_____	Pass / Fail
3 - Verify that the HIGH logic level voltage at J18-4 is 5 ± 0.5 VDC.	4.5 to 5.5 VDC	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 42 of 100	REV F
---	---	---	-----------------

5.30 Buzzer Board Interface – Inhibit Key

Test Number TCHWSCP29; Referenced Requirement HWSCP50.

5.30.1 Test Steps

- 1 - Connect an oscilloscope to J18-7 (INHIBIT_STOP_INV).
- 2 – Start ventilation with the ventilator
- 3 - Press the “Inhibit” key.
- 4 - Verify that J18-7 goes to a LOW logic level when the “Inhibit” key is pressed.
- 5 – Ensure ventilator is in ventilating
- 6 - Switch OFF the device.
- 7 - Press the “Inhibit” key.
- 8 - Verify that J18-7 goes to a LOW logic level when the “Inhibit” key is pressed.

5.30.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that J18-7 goes to a LOW logic level when the “Inhibit” key is pressed.	LOW logic level	_____	Pass / Fail
8 - Verify that J18-7 goes to a LOW logic level when the “Inhibit” key is pressed.	LOW logic level	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 43 of 100	REV F
---	---	---	-----------------

5.31 Buzzer Board Interface – Voltage Monitor

Test Number TCHWSCP B30; Referenced Requirement HWSCP B51.

5.31.1 Test Steps

- 1 - Hold down the “Inhibit” key while turning the device on to enter Maintenance mode.
- 2 - Select “Maintenance”, then “Measurement Check” from the menus.
- 3 - Use a variable power supply to apply 0 VDC to J18-5 (TESTBUZ).
- 4 - Verify that the displayed Buzzer voltage is 0 ± 0.1 VDC.
- 5 - Adjust the power supply to apply 2.5 VDC to J18-5.
- 6 - Verify that the displayed Buzzer voltage is 2.5 ± 0.1 VDC.
- 7 - Adjust the power supply to apply 5 VDC to J18-5.
- 8 - Verify that the displayed Buzzer voltage is 5 ± 0.1 VDC.

5.31.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that the displayed Buzzer voltage is 0 ± 0.1 VDC.	-0.1 to 0.1 VDC	_____	Pass / Fail
6 - Verify that the displayed Buzzer voltage is 2.5 ± 0.1 VDC.	2.4 to 2.6 VDC	_____	Pass / Fail
8 - Verify that the displayed Buzzer voltage is 5 ± 0.1 VDC.	4.9 to 5.1 VDC	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 44 of 100	REV F
---	---	---	-----------------

5.32 Remote Alarm Interface – N.C. Relay Contacts

Test Number TCHWSCP31; Referenced Requirement HWSCP52 and HWSCP54.

5.32.1 Test Steps

- 1 - Generate a LOW logic level on pin P2.14 (IC4-63) by setting the appropriate registers.
- 2 - Place a 24 VDC regulated power supply, in series with a 240 Ohm resistor and an ammeter, between J2-1 and J2-2.
- 3 - Verify that the measured current is 100 mA \pm 10%.
- 4 - Generate a HIGH logic level on pin P2.14 (IC4-63) by setting the appropriate registers.
- 5 - Verify that the measured current is 0 A.

5.32.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that the measured current is 100 mA \pm 10%.	90 mA to 110 mA	_____	Pass / Fail
5 - Verify that the measured current is 0 A.	0 A	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 45 of 100	REV F
---	---	---	-----------------

5.33 Remote Alarm Interface – N.O. Relay Contacts

Test Number TCHWSCP532; Referenced Requirement HWSCP53 and HWSCP54.

5.33.1 Test Steps

- 1 - Generate a LOW logic level on pin P2.14 (IC4-63) by setting the appropriate registers.
- 2 - Place a 24 VDC regulated power supply, in series with a 240 Ohm resistor and an ammeter, between J2-1 and J2-3.
- 3 - Verify that the measured current is 0 A.
- 4 - Generate a HIGH logic level on pin P2.14 (IC4-63) by setting the appropriate registers.
- 5 - Verify that the measured current is 100 mA \pm 10%.

5.33.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that the measured current is 0 A.	0 A	_____	Pass / Fail
5 - Verify that the measured current is 100 mA \pm 10%.	90 mA to 110 mA	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

5.34 Keyboard Interface – “Down” Key

Test Number TCHWSCP33; Referenced Requirement HWSCP55.

5.34.1 Test Steps

- 1 - Generate a LOW logic level on pin P3.0 (IC4-65) by setting the appropriate registers.
- 2 - Read the logic level on pin P3.3 (IC4-68).
- 3 - Verify that a HIGH logic level is present on P3.3 (IC4-68).
- 4 - Press the “Down” key.
- 5 - Verify that a LOW logic level is present on pin P3.3.

5.34.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that a HIGH logic level is present on P3.3 (IC4-68).	HIGH logic level	_____	Pass / Fail
5 - Verify that a LOW logic level is present on pin P3.3.	LOW logic level	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 47 of 100	REV F
---	---	---	-----------------

5.35 Keyboard Interface – “Up” Key

Test Number TCHWSCP34; Referenced Requirement HWSCP56.

5.35.1 Test Steps

- 1 - Generate a LOW logic level on pin P3.2 (IC4-67) by setting the appropriate registers.
- 2 - Read the logic level on pin P3.4 (IC4-69).
- 3 - Verify that a HIGH logic level is present on P3.4 (IC4-69).
- 4 - Press the “Up” key.
- 5 - Verify that a LOW logic level is present on pin P3.4.

5.35.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that a HIGH logic level is present on P3.4 (IC4-69).	HIGH logic level	_____	Pass / Fail
5 - Verify that a LOW logic level is present on pin P3.4.	LOW logic level	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 48 of 100	REV F
---	---	---	-----------------

5.36 Keyboard Interface – “Ventil” Key

Test Number TCHWSCP35; Referenced Requirement HWSCP57.

5.36.1 Test Steps

- 1 - Generate a LOW logic level on pin P3.1 (IC4-66) by setting the appropriate registers.
- 2 - Read the logic level on pin P3.4 (IC4-69).
- 3 - Verify that a HIGH logic level is present on pin P3.4.
- 4 - Press the “Ventil” key.
- 5 - Verify that a LOW logic level is present on pin P3.4.

5.36.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that a HIGH logic level is present on pin P3.4.	HIGH logic level	_____	Pass / Fail
5 - Verify that a LOW logic level is present on pin P3.4.	LOW logic level	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 49 of 100	REV F
---	---	---	-----------------

5.37 Keyboard Interface – “Navig” Key

Test Number TCHWSCP36; Referenced Requirement HWSCP58.

5.37.1 Test Steps

- 1 - Generate a LOW logic level on pin P3.2 (IC4-67) by setting the appropriate registers.
- 2 - Read the logic level on pin P3.3 (IC4-68).
- 3 - Verify that a HIGH logic level is present on pin P3.3.
- 4 - Press the “Navig” key.
- 5 - Verify that a LOW logic level is present on pin P3.3.

5.37.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that a HIGH logic level is present on pin P3.3.	HIGH logic level	_____	Pass / Fail
5 - Verify that a LOW logic level is present on pin P3.3.	LOW logic level	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

5.38 Keyboard Interface – “Valid” Key

Test Number TCHWSCP37; Referenced Requirement HWSCP59.

5.38.1 Test Steps

- 1 - Generate a LOW logic level on pin P3.1 (IC4-66) by setting the appropriate registers.
- 2 - Read the logic level on pin P3.3 (IC4-68).
- 3 - Verify that a HIGH logic level is present on pin P3.3.
- 4 - Press the “Valid” key.
- 5 - Verify that a LOW logic level is present on pin P3.3.

5.38.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that a HIGH logic level is present on pin P3.3.	HIGH logic level	_____	Pass / Fail
5 - Verify that a LOW logic level is present on pin P3.3.	LOW logic level	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 51 of 100	REV F
---	---	---	-----------------

5.39 Keyboard Interface – “Inhib” Key

Test Number TCHWSCP B38; Referenced Requirement HWSCP B60.

5.39.1 Test Steps

- 1 - Read the logic level on pin P3.6 (IC4-73).
- 2 - Verify that a HIGH logic level is present on pin P3.6.
- 3 - Press the “Inhib Alarm” key.
- 4 - Verify that a LOW logic level is present on pin P3.6.

5.39.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - Verify that a HIGH logic level is present on pin P3.6.	HIGH logic level	_____	Pass / Fail
4 - Verify that a LOW logic level is present on pin P3.6.	LOW logic level	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 52 of 100	REV F
---	---	---	-----------------

5.40 Keyboard Interface – Battery LED

Test Number TCHWSCP40; Referenced Requirement HWSCP39.

5.40.1 Test Steps

- 1 - Use a regulated power supply to apply 0 VDC to J7-17 (LED-BAT).
- 2 - Verify that the “Battery” LED is OFF.
- 3 - Use a regulated power supply to apply 5 VDC to J7-17 (LED-BAT).
- 4 - Verify that the “Battery” LED is ON.
- 5 - Turn the main power switch to the OFF position.
- 6 - Verify that the “Battery” LED is ON.
- 7 - Turn the main power switch to the ON position.

5.40.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - Verify that the “Battery” LED is OFF.	LED is OFF	_____	Pass / Fail
4 - Verify that the “Battery” LED is ON.	LED is ON	_____	Pass / Fail
6 - Verify that the “Battery” LED is ON.	LED is ON	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 53 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

This document may not be copied, disclosed, or used in whole or in part without the consent of COVIDIEN.

© 2008 COVIDIEN All Rights Reserved.

5.41 Keyboard Interface – DC Presence LED

Test Number TCHWSCP41; Referenced Requirement HWSCP38.

5.41.1 Test Steps

- 1 - Use a regulated power supply to apply 0 VDC to J7-16 (LED-PRESENCE-DC).
- 2 - Verify that the “DC Presence” LED is OFF.
- 3 - Use a regulated power supply to apply 5 VDC to J7-16 (LED-PRESENCE-DC).
- 4 - Verify that the “DC Presence” LED is ON.
- 5 - Turn the main power switch to the OFF position.
- 6 - Verify that the “DC Presence” LED is ON.
- 7 - Turn the main power switch to the ON position.

5.41.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - Verify that the “DC Presence” LED is OFF.	LED is OFF	_____	Pass / Fail
4 - Verify that the “DC Presence” LED is ON.	LED is ON	_____	Pass / Fail
6 - Verify that the “DC Presence” LED is ON.	LED is ON	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 54 of 100	REV F
---	---	---	-----------------

5.42 Keyboard Interface – AC Presence LED

Test Number TCHWSCP42; Referenced Requirement HWSCP37.

5.42.1 Test Steps

- 1 - Use a regulated power supply to apply 0 VDC to J7-15 (LED-PRESENCE-AC).
- 2 - Verify that the “AC Presence” LED is OFF.
- 3 - Use a regulated power supply to apply 5 VDC to J7-15 (LED-PRESENCE-AC).
- 4 - Verify that the “AC Presence” LED is ON.
- 5 - Turn the main power switch to the OFF position.
- 6 - Verify that the “AC Presence” LED is ON.
- 7 - Turn the main power switch to the ON position.

5.42.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - Verify that the “AC Presence” LED is OFF.	LED is OFF	_____	Pass / Fail
4 - Verify that the “AC Presence” LED is ON.	LED is ON	_____	Pass / Fail
6 - Verify that the “AC Presence” LED is ON.	LED is ON	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 55 of 100	REV F
---	---	---	-----------------

5.43 Keyboard Interface – Other LEDs

Test Number TCHWSCP43; Referenced Requirement HWSCP63, HWSCP64 and HWSCP65.

5.43.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start the ventilator.
- 3 - Verify that the "Ventil" LED lights at startup.
- 4 - Verify that the Red Alarm LED lights at startup.
- 5 - Verify that the Orange Alarm LED lights at startup.

5.43.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that the "Ventil" LED lights at startup.	LED is ON @ startup	_____	Pass / Fail
4 - Verify that the Red Alarm LED lights at startup.	LED is ON @ startup	_____	Pass / Fail
5 - Verify that the Orange Alarm LED lights at startup.	LED is ON @ startup	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 56 of 100	REV F
---	---	---	-----------------

5.44 Display Interface

Test Number TCHWSCP44; Referenced Requirement HWSCP66, HWSCP67 and HWSCP68.

5.44.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start the ventilator.
- 3 - Verify proper function of Display Interface.

5.44.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify proper function of Display Interface.	Displays Test Menu	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

5.45 O₂ Flow Sensor Interface (Future option)

Test Number TCHWSCP46; Referenced Requirement HWSCP71.

5.45.1 Test Steps

- 1 - Use a regulated power supply to apply 0 VDC between J20-1 and J20-3.
- 2 - Read the digital value of pin P5.10 (IC4-39, Q-O2-MEAS) by setting the appropriate registers.
- 3 - Verify that the digital value of pin P5.10 is 0 ± 11 .
- 4 - Use a regulated power supply to apply 2.5 VDC between J20-1 and J20-3.
- 5 - Read the digital value of pin P5.10 by setting the appropriate registers.
- 6 - Verify that the digital value of pin P5.10 is 512 ± 11 .
- 7 - Use a regulated power supply to apply 5 VDC between J20-1 and J20-3.
- 8 - Read the digital value of pin P5.10 by setting the appropriate registers.
- 9 - Verify that the digital value of pin P5.10 is 1023 ± 11 .

5.45.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that the digital value of pin P5.10 is 0 ± 11 .	-11 to 11	_____	Pass / Fail
6 - Verify that the digital value of pin P5.10 is 512 ± 11 .	501 to 523	_____	Pass / Fail
9 - Verify that the digital value of pin P5.10 is 1023 ± 11 .	1012 to 1034	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 58 of 100	REV F
---	---	---	-----------------

5.46 O₂ Flow Sensor Interface (Future option)

Test Number TCHWSCP47; Referenced Requirement HWSCP72.

5.46.1 Test Steps

- 1 - Connect an air compressor in series with an air flow regulator, a flow meter, and the O₂ flow sensor in the device.
- 2 - Read the digital value of pin P5.10 (IC4-39, Q-O2-MEAS) by setting the appropriate registers.
- 3 - Regulate air flow in order to measure values with the calibrated flow meter from 0 to 50 sccm, in accordance with the table. Compare the flow readings of the device to those indicated by the calibrated flow meter.
- 4 - Verify that the device readings agree with those measured by the flow meter, to within the tolerances specified in the table below:

Flow, sccm	Tolerance %
50	5
40	5
20	8
5	10
0	N/A

5.46.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - 50 ± 5% sccm	47.5 to 52.5 sccm	_____	Pass / Fail
4 - 40 ± 5% sccm	38 to 42 sccm	_____	Pass / Fail
4 - 20 ± 8% sccm	18.4 to 21.6 sccm	_____	Pass / Fail
4 - 5 ± 10% sccm	4.5 to 5.5 sccm	_____	Pass / Fail
4 - 0 ± N/A% sccm	0 sccm	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

5.47 O₂ Valve Interface – PWM and Current Drive

Test Number TCHWSCP48; Referenced Requirement HWSCP73, HWSCP74, HWSCP75 and HWSCP76.

5.47.1 Test Steps

- 1 - Generate a 14.975 kHz PWM signal with a pulsewidth of 33.38 uS on pin P7.2 by setting the appropriate registers.
- 2 - Connect an oscilloscope to pin 7.2, and measure the frequency.
- 3 - Verify that the frequency is 14.975 kHz \pm 10%.
- 4 - Generate a 14.975 kHz PWM signal with a pulsewidth of approximately 50 nS on pin P7.2 by setting the appropriate registers.
- 5 - Measure pulsewidth with the oscilloscope.
- 6 - Generate a 14.975 kHz PWM signal with a pulsewidth of approximately 75 nS on pin P7.2 by setting the appropriate registers.
- 7 - Measure pulsewidth with the oscilloscope.
- 8 - Verify that the difference between first and second pulsewidth measurement is 25 \pm 5 nS.
- 9 - Unplug J21 (the O₂ Valve connector).
- 10 - Generate a HIGH logic level signal on P7.2 (IC4-21, CD-O2) by setting the appropriate microcontroller registers.
- 11 - Connect a 110 Ohm, 5 Watt resistor, in series with an ammeter, between J21-1 and J21-2.
- 12 - Verify that the current through the resistor \geq 190 mA.

5.47.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that the frequency is 14.975 kHz \pm 10%.	13.4775 to 16.4725 kHz	_____	Pass / Fail
5 - Measure pulsewidth with the oscilloscope.	N/A	_____	N/A
7 - Measure pulsewidth with the oscilloscope.	N/A	_____	N/A
8 - Verify that the difference between first and second pulsewidth measurement is 25 \pm 5 nS.	20 to 30 nS	_____	Pass / Fail
12 - Verify that the current through the resistor is \geq 190 mA.	\geq 190 mA	_____	Pass / Fail

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 60 of 100	REV F
---	---	---	-----------------

Comments: _____

Test Operator Signature/Date: _____

5.48 Security Back-up Buzzer Interface – Secondary Alarm

Test Number TCHWSCP49; Referenced Requirement HWSCP77.

5.48.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Run "Buzzer Test", and skip to the "Security Buzzer" test..
- 3 - Run the "Security Buzzer Test".
- 4 - Verify that the secondary audible alarm sounds when the device resets.
- 5 - Verify that the secondary audible alarm's Sound Pressure Level is at least 80 dB (A-weighted) at a distance of 10 cm.

5.48.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that the secondary audible alarm sounds when the device resets.	Alarm sounds	_____	Pass / Fail
5 - Verify that the secondary audible alarm's Sound Pressure Level is at least 80 dB (A-weighted) at a distance of 10 cm.	80 dB	_____	N/A

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 62 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

5.49 Security Back-up Buzzer Interface – Primary Alarm

Test Number TCHWSCP50; Referenced Requirement HWSCP78.

5.49.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start ventilator in "Partial Tests" mode.
- 3 - Run the "Buzzer Test".
- 4 - Verify that the primary audible alarm sounds when the device resets.
- 5 - Verify that the primary audible alarm's Sound Pressure Level is at least 90 dB (A-weighted) at a distance of 10 cm.

5.49.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that the primary audible alarm sounds when the device resets.	Alarm sounds	_____	Pass / Fail
5 - Verify that the primary audible alarm's Sound Pressure Level is at least 90 dB (A-weighted) at a distance of 10 cm.	90 dB	_____	N/A

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 63 of 100	REV F
---	---	---	-----------------

5.50 RAM Memory

Test Number TCHWSCP51; Referenced Requirement HWSCP79.

5.50.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start the ventilator in "Partial Tests" mode.
- 3 - Run "RAM Memory" test.
- 4 - Verify that the result of the RAM Memory test is "OK".

5.50.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that the result of the RAM Memory test is "OK".	OK	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 64 of 100	REV F
---	---	---	-----------------

5.51 Event Memory

Test Number TCHWSCP52; Referenced Requirement HWSCP80.

5.51.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start the ventilator in "Partial Tests" mode.
- 3 - Run "Event Memory" test.
- 4 - Verify that the result of the Events Memory test is "OK".

5.51.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that the result of the Events Memory test is "OK".	OK	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 65 of 100	REV F
---	---	---	-----------------

5.52 Event Memory – Erase/Program Cycles

Test Number TCHWSCP53; Referenced Requirement HWSCP81.

5.52.1 Test Steps

- 1 - Refer to the specifications for the non-volatile memory.
- 2 - By inspection of the specifications, verify that the memory will support at least 1,00,000 erase/program cycles.

5.52.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - By inspection verify that the memory will support at least 1,000,000 erase and program cycles.	>= 1,00,000	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 66 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

5.53 Ventilation Settings Memory

Test Number TCHWSCP82; Referenced Requirement HWSCP82.

5.53.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start the ventilator in "Partial Tests" mode.
- 3 - Run "Parameters EEPROM Memory" test.
- 4 - Verify that the result of the test is "OK".

5.53.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that the result of the Parameters EEPROM Memory test is "OK".	OK	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 67 of 100	REV F
---	---	---	-----------------

5.54 Ventilation Settings Memory – Erase/Program Cycles

Test Number TCHWSCP83; Referenced Requirement HWSCP83.

5.54.1 Test Steps

- 1 - Refer to the specifications for the non-volatile memory.
- 2 - By inspection of the specifications, verify that the memory will support at least 1,000,000 erase/program cycles.

5.54.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - By inspection verify that the memory will support at least 1,000,000 erase and program cycles.	>= 1,000,000	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 68 of 100	REV F
---	---	---	-----------------

5.55 Monitoring Memory

Test Number TCHWSCP56; Referenced Requirement HWSCP84.

5.55.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start the ventilator in "Partial Tests" mode.
- 3 - Run "Monitoring Memory" test.
- 4 - Verify that the result of the Monitoring Memory test is "OK".

5.55.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that the result of the Monitoring Memory test is "OK".	OK	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 69 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

5.56 Monitoring Memory – Erase/Program Cycles

Test Number TCHWSCP857; Referenced Requirement HWSCP85.

5.56.1 Test Steps

- 1 - Refer to the specifications for the non-volatile memory.
- 2 - By inspection of the specifications, verify that the memory will support at least 100,000 erase/program cycles.

5.56.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - By inspection verify that the memory will support at least 100,000 erase and program cycles.	>= 100,000	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____


COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 70 of 100	REV F
---	---	---	-----------------

5.57 Watchdog – Power Supervisor

Test Number TCHWSCP B58; Referenced Requirements HWSCP B89, HWSCP B90, HWSCP B91, and HWSCP B92.

5.57.1 Test Steps


- 1 - Turn off the unit and disconnect it from AC power.
- 2 - Disconnect one end of the cable between the CPU Board and the Power Management Board.
- 3 - Connect a variable, programmable power supply to the input of IC26 (the pin closest to the long edge of the CPU board) and Ground, to substitute for the +5V supply voltage.
- 4 - Connect oscilloscope channel 1 to T3-Drain (RESET). The Drain lead is on the side of T3 that has only one lead.
- 5 - Slowly reduce the power supply voltage until the RESET signal on oscilloscope channel 1 goes to a HIGH logic level, and record the voltage at which this occurs.
- 6 - Verify that the voltage reading is between 4.5 and 4.75 VDC.
- 7 - Slowly increase the power supply voltage until the RESET signal on oscilloscope channel 1 goes LOW again, and record the voltage at which this occurs.
- 8 - Verify that the voltage reading in step 7 is at least 0.01 VDC higher than the reading from step 6.
- 9 - Adjust the power supply to produce a repeating voltage waveform that is at 5 VDC for 4.9 seconds, and at 4 VDC for 0.1 seconds.
- 10 - Connect oscilloscope channel 1 to the power supply output, and connect channel 2 to T3-Drain (RESET).
- 11 - Verify that the pulsewidth of the 5 V portion of the RESET signal on scope channel 2 has a pulsewidth of at least 140 ms, and no more than 280 ms, as shown in the diagram below.



- 12 - While observing the oscilloscope, increase the length of time that the power supply voltage is at 4 V to 400 ms.

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 71 of 100	REV F
---	---	---	-----------------

13 - Verify that RESET continues to stay high for 140 to 280 ms after Vcc returns to normal, as shown in the diagram below.



14 - Turn off the power supply and disconnect it from IC12-2, and reconnect the cable between the CPU Board and Power Management Board.

5.57.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
6 - Verify that the voltage at which RESET goes high is between 4.5 and 4.75 VDC.	4.5 to 4.75 VDC	_____	Pass / Fail
8 - Verify that the voltmeter reading is at least 0.01 VDC higher than the previous reading.	≥ 0.01 VDC higher than previous reading	_____	Pass / Fail
11 - Verify that the 5 V portion of the RESET pulse is between 140 and 280 ms in duration.	140 to 280 ms	_____	Pass / Fail
13 - Verify that RESET goes low again between 140 and 280 ms after Vcc returns to normal.	140 to 280 ms	_____	Pass / Fail

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 72 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

This document may not be copied, disclosed, or used in whole or in part without the consent of COVIDIEN.

© 2008 COVIDIEN All Rights Reserved.

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 73 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

This document may not be copied, disclosed, or used in whole or in part without the consent of COVIDIEN.

© 2008 COVIDIEN All Rights Reserved.

5.58 Watchdog – Timeout interval

Test Number TCHWSCP859; Referenced Requirement HWSCP86 and HWSCP88.

5.58.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Connect oscilloscope channel 1 to IC12-6 (WDOG), and channel 2 to T3-Drain (RESET).
- 3 - Start ventilator in "Partial Tests" mode.
- 4 - Run "Watchdog Function" test, and capture the signals displayed on the oscilloscope.
- 5 - Verify that the delay between the last level transition of oscilloscope channel 1 (WDOG), and oscilloscope channel 2 (RESET) going high is between 1.0 and 2.25 seconds.
- 6 - Verify that the pulsewidth of oscilloscope ch. 2 (RESET) is between 140 and 280 mS.

5.58.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
5 - Verify that the delay between the last level transition of oscilloscope ch. 1 (WDOG), and oscilloscope ch. 2 (RESET) going high is between 1.0 and 2.25 seconds.	1.0 to 2.25 seconds	_____	Pass / Fail
6 - Verify that the pulsewidth of oscilloscope ch. 2 (RESET) is between 140 and 280 mS.	140 to 280 mS	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 74 of 100	REV F
---	---	---	-----------------

5.59 Watchdog – PFI Trap Function

Test Number TCHWSCP B73; Referenced Requirement HWSCP B93.

5.59.1 Test Steps

- 1 - Disconnect the cable from the Power Management Board to the CPU board at the Power Management Board end.
- 2 - Connect the positive lead of a variable power supply, adjusted to 24 VDC, to J7-3 (+24VUTIL). Connect the negative power supply lead J7-24 (Ground). Monitor the power supply's voltage with its built-in meter, or with an external voltmeter.
- 3 - Wait until the unit proceeds to the ventilation status screen. Disregard any error messages that may appear in the status area.
- 4 - Reduce the power supply voltage to 14 VDC. Do not allow the voltage to go below 14 VDC during this adjustment.
- 5 - While observing the unit's LCD display, slowly decrease the power supply voltage until the LCD display goes black.
- 6 - Verify that the power supply voltage is between 12.19 VDC and 13.83 VDC.
- 7 - Allow the unit to reset, to silence the power failure alarm.

5.59.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
6 - Verify that the power supply voltage measures from 12.19 to 13.83 VDC.	12.19 to 13.83 VDC	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 75 of 100	REV F
---	---	---	-----------------

5.60 Clock – Functional Test

Test Number TCHWSCP60; Referenced Requirement HWSCP694.

5.60.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start ventilator in "Partial Tests" mode.
- 3 - Run the "Real Time Clock" test.
- 4 - Verify that the result of the Real Time Clock test is "OK".

5.60.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
4 - Verify that the result of the Real Time Clock function test is "OK".	OK	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 76 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

This document may not be copied, disclosed, or used in whole or in part without the consent of COVIDIEN.

© 2008 COVIDIEN All Rights Reserved.

5.61 Clock – Battery Timekeeping Mode

Test Number TCHWSCP61; Referenced Requirement HWSCP65.

5.61.1 Test Steps

- 1 - Refer to the specifications for the button cell battery and DS1305E clock.
- 2 - By inspection verify that the button cell's capacity is sufficient to allow the DS1305E clock chip to operate in battery timekeeping mode for a minimum of three years, even if the device is not supplied with power.

5.61.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - By datasheet inspection, verify that the button cell's capacity is sufficient to allow the DS1305E clock chip to operate in battery timekeeping mode for a minimum of three years, even if the device is not supplied with power.	> 3 years	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 77 of 100	REV F
---	---	---	-----------------

5.62 Device Model ID

Test Number TCHWSCP62; Referenced Requirement HWSCP696.

5.62.1 Test Steps

- 1 - Read the logic level of pin P8.2 (IC4-11).
- 2 - Verify that a HIGH logic level is present on pin P8.2.
- 3 - Read the logic level of pin P8.3 (IC4-12).
- 4 - Verify that a HIGH logic level is present on pin P8.3.

5.62.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail / N/A
2 - Verify that a HIGH logic level is present on pin P8.2.	HIGH	_____	Pass / Fail
4 - Verify that a HIGH logic level is present on pin P8.3.	HIGH	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 78 of 100	REV F
---	---	---	-----------------

5.63 Inspiratory Flow Measurement

Test Number TCHWSCP63; Referenced Requirement HWSCP97 and HWSCP98.

5.63.1 Test Steps

- 1 - Connect an air compressor in series with an air flow regulator, a flow meter and the inhalation flow sensor.
- 2 - Read the digital value of the signal on pin P5.1 (IC4-28, QI-MEAS) by setting the appropriate registers.
- 3 - Regulate air flow in order to measure values with the flow meter from 0 to 1000 sccm, in increments of 100 sccm.
- 4 - Verify that the digital value matches following table values:

Flow sccm	Nominal Voltage V (from sensor data-sheet)	Tolerance \pm V DC	Tolerance %. Not Calibrated. Sensor only.	Result
1000	5.00	0.5000	10.00%	
900	4.90	0.5050	10.31%	
800	4.80	0.5100	10.63%	
700	4.66	0.5130	11.01%	
600	4.42	0.5110	11.56%	
500	4.18	0.5090	12.18%	
400	3.82	0.5010	13.12%	
300	3.41	0.4605	13.50%	
200	2.96	0.4180	14.12%	
100	2.30	0.3550	15.43%	
0	1.00	0.2500	25.00%	

Flow l/min	Tolerance % Calibrated. With whole supply chain measure.	Result
145	5.00%	
135	5.00%	
90	5.00%	
60	5.00%	
37	5.00%	
12	8.00%	
5	10.00%	
0	N/A	

5.63.2 Results

Description Of Expected Result	Expected Result	Pass / Fail
Verify inspiratory air flow sensor's accuracy.	Per above tables.	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

5.64 Exhalation Flow Measurement

Test Number TCHWSCP64; Referenced Requirement HWSCP69 and HWSCP100.

5.64.1 Test Steps

- 1 - Connect an air compressor in series with an air flow regulator, a flow meter and the exhalation flow sensor.
- 2 - Read the digital value of the signal on pin P5.7 (IC4-28, QI-MEAS) by setting the appropriate registers.
- 3 - Regulate air flow in order to measure values with the flow meter from 0 to 1000 sccm, in increments of 100 sccm.
- 4 - Verify that the digital value matches following table values:

Flow sccm	Nominal Voltage V (from sensor data-sheet)	Tolerance \pm V DC	Tolerance %. Not Calibrated. Sensor only.	Result
1000	5.00	0.5000	10.00%	
900	4.90	0.5050	10.31%	
800	4.80	0.5100	10.63%	
700	4.66	0.5130	11.01%	
600	4.42	0.5110	11.56%	
500	4.18	0.5090	12.18%	
400	3.82	0.5010	13.12%	
300	3.41	0.4605	13.50%	
200	2.96	0.4180	14.12%	
100	2.30	0.3550	15.43%	
0	1.00	0.2500	25.00%	

Flow l/min	Tolerance % Calibrated. With whole supply chain measure.	Result
145	5.00%	
135	5.00%	
90	5.00%	
60	5.00%	
37	5.00%	
12	8.00%	
5	10.00%	
0	N/A	

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 81 of 100	REV F
---	---	---	-----------------

5.64.2 Results

Description Of Expected Result	Expected Result	Pass / Fail
Verify expiratory air flow sensor's accuracy.	Per above tables.	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 82 of 100	REV F
---	---	---	-----------------

5.65 Proximal Pressure Measurement

Test Number TCHWSCP65; Referenced Requirement HWSCP101.

5.65.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start ventilator in "Partial Tests" mode.
- 3 - Run "Pressure Sensor" test.
- 4 - Use a pressure meter to measure the pressure applied to the proximal pressure sensor.
- 5 - Apply a -10 mbar constant proximal pressure (suction) to the proximal pressure sensor with a syringe.
- 6 - Verify that the measured pressure matches the proximal pressure displayed by the device, to within 1.5 mbar.
- 7 - Apply a 50 mbar constant proximal pressure with a syringe.
- 8 - Verify that the measured pressure matches the proximal pressure displayed by the device, to within 1.5 mbar.
- 9 - Apply a 100 mbar constant proximal pressure with a syringe.
- 10 - Verify that the measured pressure matches the proximal pressure displayed by the device, to within 1.5 mbar.

5.65.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
6 - Verify that the measured pressure matches the proximal pressure displayed by the device, to within 1.5 mbar.	-11.5 to -8.5 mbar	_____	Pass / Fail
8 - Verify that the measured pressure matches the proximal pressure displayed by the device, to within 1.5 mbar.	48.5 to 51.5 mbar	_____	Pass / Fail
10 - Verify that the measured pressure matches the proximal pressure displayed by the device, to within 1.5 mbar.	98.5 to 101.5 mbar	_____	Pass / Fail

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 83 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

Comments: _____

Test Operator Signature/Date: _____

5.66 Internal Pressure Measurement

Test Number TCHWSCP66; Referenced Requirement HWSCP102.

5.66.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start ventilator in "Partial Tests" mode.
- 3 - Run "Pressure Sensor" test.
- 4 - Use a pressure meter to measure the pressure applied to the inspiratory pressure sensor.
- 5 - Apply a -10 mbar constant internal pressure (suction) to the inspiratory pressure sensor with a syringe.
- 6 - Verify that the measured internal pressure matches that displayed by the device, to within 1.5 mbar
- 7 - Apply a 50 mbar constant internal pressure with a syringe.
- 8 - Verify that the measured internal pressure matches that displayed by the device, to within 1.5 mbar.
- 9 - Apply a 100 mbar constant internal pressure with a syringe.
- 10 - Verify that the measured internal pressure matches that displayed by the device, to within 1.5 mbar.

5.66.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
6 - Verify that the measured pressure matches the internal pressure displayed by the device, to within 1.5 mbar.	-11.5 to -8.5 mbar	_____	Pass / Fail
8 - Verify that the measured pressure matches the internal pressure displayed by the device, to within 1.5 mbar.	48.5 to 51.5 mbar	_____	Pass / Fail
10 - Verify that the measured pressure matches the internal pressure displayed by the device, to within 1.5 mbar.	98.5 to 101.5 mbar	_____	Pass / Fail

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 85 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

Comments: _____

Test Operator Signature/Date: _____

5.67 Expiratory Valve Pressure Measurement

Test Number TCHWSCP67; Referenced Requirement HWSCP103.

5.67.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start ventilator in "Partial Tests" mode.
- 3 - Run "Pressure Sensor" test.
- 4 - Use a pressure meter to measure the pressure applied to the expiratory valve pressure sensor.
- 5 - Apply a -10 mbar constant internal pressure (suction) to the expiratory valve pressure sensor with a syringe.
- 6 - Verify that the measured internal pressure matches that displayed by the device, to within 1.5 mbar
- 7 - Apply a 50 mbar constant pressure with a syringe.
- 8 - Verify that the measured internal pressure matches that displayed by the device, to within 1.5 mbar.
- 9 - Apply a 100 mbar constant pressure with a syringe.
- 10 - Verify that the measured expiratory valve pressure matches that displayed by the device, to within 1.5 mbar.

5.67.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
6 - Verify that the measured pressure matches the proximal pressure displayed by the device, to within 1.5 mbar.	-11.5 to -8.5 mbar	_____	Pass / Fail
8 - Verify that the measured pressure matches the proximal pressure displayed by the device, to within 1.5 mbar.	48.5 to 51.5 mbar	_____	Pass / Fail
10 - Verify that the measured pressure matches the proximal pressure displayed by the device, to within 1.5 mbar.	98.5 to 101.5 mbar	_____	Pass / Fail

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 87 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

Comments: _____

Test Operator Signature/Date: _____

5.68 Atmospheric Pressure Measurement

Test Number TCHWSCP68; Referenced Requirement HWSCP104.

5.68.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start ventilator in "Partial Tests" mode.
- 3 - Run "Pressure Sensor" test.
- 4 - Use an absolute pressure meter to measure absolute pressure.
- 5 - Verify that the measured absolute pressure matches device displayed ambient pressure value within 15 mbar.
- 6 - Apply a negative pressure on absolute pressure sensor with a syringe until the device displays an ambient pressure of 600 mbar.
- 7 - Use an absolute pressure meter to measure this pressure.
- 8 - Verify that the measured pressure is equal to 600 ± 15 mbar.
- 9 - Apply a positive pressure on absolute pressure sensor with a syringe until the device displays an ambient pressure of 1100 mbar.
- 10 - Use an absolute pressure meter to measure this pressure.
- 11 - Verify that the measured pressure is equal to 1100 ± 15 mbar.

5.68.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
5 - Verify that the measured absolute pressure matches device displayed ambient pressure value within 10 mbar.	(Ambient – 15) to (Ambient + 15) mbar	_____	Pass / Fail
8 - Verify that the measured pressure is equal to 600 ± 15 mbar.	585 to 615 mbar	_____	Pass / Fail
11 - Verify that the measured pressure is equal to 1100 ± 15 mbar.	1085 to 1115 mbar	_____	Pass / Fail

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 89 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

Comments: _____

Test Operator Signature/Date: _____

5.69 O₂ Pressure Measurement (Future option)

Test Number TCHWSCP69; Referenced Requirement HWSCP105.

5.69.1 Test Steps

- 1 - Download approved "CPU Test Software" version to the CPU Board, and approved "Soft carte alim pour test CPU" software version to the Power Supply Management board.
- 2 - Start ventilator in "Partial Tests" mode.
- 3 - Run "Pressure Sensor" test.
- 4 - Use a pressure meter to measure O₂ pressure.
- 5 - Apply a 0 bar constant O₂ pressure with a syringe.
- 6 - Verify that the measured O₂ pressure matches device displayed O₂ pressure value to within 150 mbar.
- 7 - Apply a 3.5 bar constant O₂ pressure with a syringe.
- 8 - Verify that the measured O₂ pressure matches device displayed O₂ pressure value to within 150 mbar.
- 9 - Apply a 6.5 bar constant O₂ pressure with a syringe.
- 10 - Verify that the measured O₂ pressure matches device displayed O₂ pressure value to within 150 mbar.

5.69.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
6 - Verify that the device-displayed O ₂ pressure matches the measured O ₂ pressure value to within 150 mbar.	(Measured – 150) to (Measured + 150) mbar	_____	Pass / Fail
8 - Verify that the measured O ₂ pressure matches device displayed O ₂ pressure value to within 150 mbar.	3.35 to 3.65 mbar	_____	Pass / Fail
10 - Verify that the measured O ₂ pressure matches device displayed O ₂ pressure value to within 150 mbar.	6.35 to 6.65 mbar	_____	Pass / Fail

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 91 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

Comments: _____

Test Operator Signature/Date: _____


COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 92 of 100	REV F
---	---	---	-----------------

5.70 PCB – Size

Test Number TCHWSCP70; Referenced Requirement HWSCP106, HWSCP107, HWSCP108, HWSCP109, HWSCP110 and HWSCP111.

5.70.1 Test Steps

1 - Measure the dimensions of the CPU Board's PCB, as shown in the following diagram:



- 2 - Verify that the length (L1) is 190 ± 1 mm.
- 3 - Verify that the length (L2) is 58.5 ± 1 mm.
- 4 - Verify that the length (L3) is 96.5 ± 1 mm.
- 5 - Verify that the width (W1) is 100 ± 1 mm.
- 6 - Verify that the width (W2) is 20 ± 1 mm.
- 7 - Verify that the thickness (T) is 1.6 ± 0.5 mm.

5.70.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - Verify that the length (L1) is 190 ± 1 mm	189 to 191 mm	_____	Pass / Fail
3 - Verify that the length (L2) is 58.5 ± 1 mm.	57.5 to 59.5 mm	_____	Pass / Fail
4 - Verify that the length (L3) is 96.5 ± 1 mm.	95.5 to 97.5 mm	_____	Pass / Fail
5 - Verify that the width (W1) is 100 ± 1 mm.	99 to 101 mm	_____	Pass / Fail
6 - Verify that the width (W2) is 20 ± 1 mm.	19 to 21 mm	_____	Pass / Fail
7 - Verify that the thickness (T) is 1.6 ± 0.5 mm.	1.1 to 2.1 mm	_____	Pass / Fail

COVIDIEN
6135 Gunbarrel Avenue
Boulder, CO 80301

TITLE:
Test Protocol, CPU Board, PB540,
PB560

DOCUMENT NUMBER

10027876

SHEET 93 of 100

REV

F

Comments: _____

Test Operator Signature/Date: _____

5.71 PCB – Mounting

Test Number TCHWSCP71; Referenced Requirement HWSCP112.

5.71.1 Test Steps

- 1 - Measure the diameter of the three PCB mounting holes, shown in the PCB Size test case.
- 2 - Verify that the holes are 3.5 ± 0.2 mm in diameter.

5.71.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - Verify that the first hole is 3.5 mm in diameter.	3.3 to 3.7 mm	_____	Pass / Fail
2 - Verify that the second hole is 3.5 mm in diameter.	3.3 to 3.7 mm	_____	Pass / Fail
2 - Verify that the third hole is 3.5 mm in diameter.	3.3 to 3.7 mm	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 95 of 100	REV F
---	---	---	-----------------

5.72 PCB – Labeling

Test Number TCHWSCP72; Referenced Requirement HWSCP114.

5.72.1 Test Steps

1 - Inspect the markings on the PCB.

2 - Verify that the PCB is labeled with its Name, Part Number, and Revision, and that the markings are human-readable.

5.72.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
2 - Verify that the PCB is labeled with its Name, and that the markings are human-readable.	Name	_____	Pass / Fail
2 - Verify that the PCB is labeled with its Part Number, and that the markings are human-readable.	Part Number	_____	Pass / Fail
2 - Verify that the PCB is labeled with its Revision, and that the markings are human-readable.	Revision	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 96 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

5.73 Printed Circuit Board Fabrication Standards

Test Number TCHWSCP73; Referenced Requirement HWSCP113.

Verify that the CPU Card is specified to be fabricated according to the following standards and exceptions:

- IPC-A-610D Class III: 2005, Acceptability for Electronic Assemblies.
- Any exceptions to these standards are listed on the fabrication drawings.

5.73.1 Test Steps

- 1 - Verify that specification or statement exists in the CPU Card's fabrication documentation set indicating that it must comply with IPC-610D Class III: 2005 and that exceptions, if any, to this standard are stated. Attach a copy of the document page(s), with the location highlighted of the specification(s) or statement(s) found providing evidence that this requirement has been met.

5.73.2 Test Results

CPU Card fabrication documentation set (in Agile):

Part Number: _____ Revision _____

Description of Expected Result	Expected Result	Actual Result	Pass / Fail
(1) CPU Card Assembly drawing states PCB meets IPC-610D Class III: 2005 and exceptions, if any, are noted.	PCB Fab Dwg includes statement that it meets IPC-610D Class III: 2005 and exceptions, if any, are noted	_____	Pass / Fail

Attach copies of pages, with corresponding statements highlighted, where the above results were found.

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 97 of 100	REV F
---	---	---	-----------------

5.74 FiO2 Measurement (PB560, Future option)

Test Number TCHWSCP B74; Referenced Requirement HWSCP B69.

5.74.1 Test Steps

- 1 – Do a short-circuit between J11-1 and J11-2 to apply 0 VDC.
- 2 - Read the digital value of pin P5.5 (IC4-32, FiO2-MEAS) by setting the appropriate registers.
- 3 - Verify that the digital value of pin P5.5 is 0 ± 10 .
- 4 - Use a regulated power supply to apply 33mVDC between J11-1 and J11-2.
- 5 - Read the digital value of pin P5.5 by setting the appropriate registers.
- 6 - Verify that the digital value of pin P5.5 is 500 ± 10 .
- 7 - Use a regulated power supply to apply 66mVDC between J11-1 and J11-2.
- 8 - Read the digital value of pin P5.5 by setting the appropriate registers.
- 9 - Verify that the digital value of pin P5.5 is 1000 ± 10 .

5.74.2 Results

Description Of Expected Result	Expected Result	Actual Result	Pass / Fail
3 - Verify that the digital value of pin P5.5 is 0 ± 10 .	-10 to 10	_____	Pass / Fail
6 - Verify that the digital value of pin P5.5 is 500 ± 10 .	490 to 510	_____	Pass / Fail
9 - Verify that the digital value of pin P5.5 is 1000 ± 10 .	990 to 1010	_____	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 98 of 100	REV F
---	---	---	-----------------

TEST SUMMARY

Test Cases	Requirements	Pass / Fail
TCHWSCP B1	HWSCP B1	Pass / Fail
TCHWSCP B2	HWSCP B2	Pass / Fail
TCHWSCP B3	HWSCP B3	Pass / Fail
TCHWSCP B4	HWSCP B4	Pass / Fail
TCHWSCP B5	HWSCP B8, HWSCP B9, HWSCP B10	Pass / Fail
TCHWSCP B6	HWSCP B5, HWSCP B6, HWSCP B7	Pass / Fail
TCHWSCP B7	HWSCP B11, HWSCP B12, HWSCP B13	Pass / Fail
TCHWSCP B8	HWSCP B14, HWSCP B15, HWSCP B16	Pass / Fail
TCHWSCP B9	HWSCP B19, HWSCP B20	Pass / Fail
TCHWSCP B10	HWSCP B17, HWSCP B18	Pass / Fail
TCHWSCP B11	HWSCP B21, HWSCP B22	Pass / Fail
TCHWSCP B12	HWSCP B23, HWSCP B25, HWSCP B26	Pass / Fail
TCHWSCP B13	HWSCP B24, HWSCP B27, HWSCP B28	Pass / Fail
TCHWSCP B14	HWSCP B29, HWSCP B30	Pass / Fail
TCHWSCP B15	HWSCP B31, HWSCP B32	Pass / Fail
TCHWSCP B16	HWSCP B33	Pass / Fail
TCHWSCP B17	HWSCP B34	Pass / Fail
TCHWSCP B18	HWSCP B35	Pass / Fail
TCHWSCP B19	HWSCP B36	Pass / Fail
TCHWSCP B20	HWSCP B43	Pass / Fail
TCHWSCP B22	HWSCP B37, HWSCP B38, HWSCP B39	Pass / Fail
TCHWSCP B23	HHWSCP B40	Pass / Fail
TCHWSCP B24	HWSCP B41	Pass / Fail
TCHWSCP B25	HWSCP B42	Pass / Fail
TCHWSCP B26	HWSCP B45	Pass / Fail
TCHWSCP B27	HWSCP B46, HWSCP B47, HWSCP B48	Pass / Fail
TCHWSCP B28	HWSCP B49	Pass / Fail
TCHWSCP B29	HWSCP B50	Pass / Fail
TCHWSCP B30	HWSCP B51	Pass / Fail
TCHWSCP B31	HWSCP B52, HWSCP B54	Pass / Fail
TCHWSCP B32	HWSCP B53, HWSCP B54	Pass / Fail
TCHWSCP B33	HWSCP B55	Pass / Fail
TCHWSCP B34	HWSCP B56	Pass / Fail
TCHWSCP B35	HWSCP B57	Pass / Fail
TCHWSCP B36	HWSCP B58	Pass / Fail
TCHWSCP B37	HWSCP B59	Pass / Fail
TCHWSCP B38	HWSCP B60	Pass / Fail
TCHWSCP B40	HWSCP B39	Pass / Fail
TCHWSCP B41	HWSCP B38	Pass / Fail
TCHWSCP B42	HWSCP B37	Pass / Fail
TCHWSCP B43	HWSCP B63, HWSCP B64, HWSCP B65	Pass / Fail
TCHWSCP B44	HWSCP B66, HWSCP B67, HWSCP B68	Pass / Fail
TCHWSCP B46	HWSCP B71	Pass / Fail
TCHWSCP B47	HWSCP B72	Pass / Fail
TCHWSCP B48	HWSCP B73, HWSCP B74, HWSCP B75, HWSCP B76	Pass / Fail
TCHWSCP B49	HWSCP B77	Pass / Fail
TCHWSCP B50	HWSCP B78	Pass / Fail
TCHWSCP B51	HWSCP B79	Pass / Fail
TCHWSCP B52	HWSCP B80	Pass / Fail

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 99 of 100	REV F
---	---	---	-----------------

PROPRIETARY & CONFIDENTIAL

FOR INTERNAL USE ONLY

TCHWSCP53	HWSCP81	Pass / Fail
TCHWSCP54	HWSCP82	Pass / Fail
TCHWSCP55	HWSCP83	Pass / Fail
TCHWSCP56	HWSCP84	Pass / Fail
TCHWSCP57	HWSCP85	Pass / Fail
TCHWSCP58	HWSCP89, HWSCP90, HWSCP91, HWSCP92	Pass / Fail
TCHWSCP59	HWSCP86, HWSCP88	Pass / Fail
TCHWSCP73	HWSCP93	Pass / Fail
TCHWSCP60	HWSCP94	Pass / Fail
TCHWSCP61	HWSCP95	Pass / Fail
TCHWSCP62	HWSCP96	Pass / Fail
TCHWSCP63	HWSCP97, HWSCP98	Pass / Fail
TCHWSCP64	HWSCP99, HWSCP100	Pass / Fail
TCHWSCP65	HWSCP101	Pass / Fail
TCHWSCP66	HWSCP102	Pass / Fail
TCHWSCP67	HWSCP103	Pass / Fail
TCHWSCP68	HWSCP104	Pass / Fail
TCHWSCP69	HWSCP105	Pass / Fail
TCHWSCP70	HWSCP106, HWSCP107, HWSCP108, HWSCP109, HWSCP110, HWSCP111	Pass / Fail
TCHWSCP71	HWSCP112	Pass / Fail
TCHWSCP72	HWSCP114	Pass / Fail
TCHWSCP73	HWSCP93, HWSCP113	Pass / Fail
TCHWSCP74	HWSCP69	Pass / Fail
TCHWSCP76	HWSCP117	Pass / Fail
TCHWSCP77	HWSCP44	Pass / Fail

Comments: _____

Test Operator Signature/Date: _____

Test Operator Printed Name: _____

Review Engineer Signature/Date: _____

Review Engineer Printed Name: _____

COVIDIEN 6135 Gunbarrel Avenue Boulder, CO 80301	TITLE: Test Protocol, CPU Board, PB540, PB560	DOCUMENT NUMBER 10027876 SHEET 100 of 100	REV F
---	---	--	-----------------